www.medela.ca
info@medela.ca
(800) 435-8316 or (905) 608-7272
Canada
Mississauga, ON, L5L 0A1,
4160 Stadelview Crescent, #8
Medela Canada Inc.

www.medela.us
customer.service@medela.com
(800) 435-8316 or (815) 363-1166
USA
McHenry, IL 60050,
1101 Corporate Drive
Medela LLC

Swing Maxi™ Hands-free

Electric breast pump
Instructions for use

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





13

1. Important safety information

! Please read all instructions before using this product. Keep these instructions for future reference.

Warnings identify all instructions that are important for safety. Failure to observe these instructions can lead to personal injury or damage to the product. The following symbols and signal words show the significance of the warnings:

! WARNING

Can lead to serious injury or death.

NOTICE

Can lead to material damage (not related to personal injury).

! CAUTION

Can lead to minor injury.

i Information

Useful or important information that is not related to safety.

When using electrical products, especially when children are present, basic safety precautions must always be followed.

! WARNING

To avoid fire, electric shock, or serious burns:

- The breast pump should never be left unattended when plugged into a power source.
- Never operate an electrical device if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. If damage is found, immediately discontinue use of power adaptor and call Medela Customer Service.
- Do not use an electrical product that has been exposed to water or other liquids including:
 - Do not use while bathing or showering.
 - Never place or drop into water or other liquids.
 - Do not run water over the breast pump.
 - If an electrical product has been exposed to water or other liquids, do not touch and if possible unplug the device from electrical outlet, turn off and contact the manufacturer.
- Do not expose the breast pump or battery to excessive temperatures.
- Only use power adaptor and cable that comes with the breast pump. An inadequate power source or cable may result in a fire hazard, electric shock or malfunction equipment due to overvoltage.

To avoid health risks and reduce the risk of severe injury:

- This is a single-user product. Powered breast pumps that are designed for single users should never be rented or shared. Use by more than one person may present a health risk and voids the warranty.

- Do not use the breast pump while operating a moving vehicle.
- Pumping may induce contractions or premature labour. Do not pump until after giving birth. If you become pregnant while breastfeeding or breast pumping, consult with a licensed healthcare professional before continuing.
- If infected with Hepatitis B, Hepatitis C, or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- Wash all parts that come into contact with your breast and breast milk after every use.
- The tubing poses a risk of strangulation when being played with and small parts may be swallowed by toddlers.
- Close supervision is necessary when the breast pump or accessories are used in the vicinity of children.
- This appliance is not intended for use by persons (including children up to 14 years) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children must not play with the appliance.
- Cleaning and maintenance must not be performed by children without supervision.

CAUTION

To avoid fire or burns:

- Make sure the voltage of the power adaptor is compatible with the power source. See **chapter 18** for technical specifications.
- Always unplug electrical products immediately after use except for when charging.
- Do not place or store the product where it can fall into water.
- The breast pump and accessories are not heat-resistant. Keep away from heated surfaces or open flames.
- Do not use near flammable materials.

To avoid health risks and reduce the risk of injury:

- Before each use visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- This product cannot be serviced or repaired. Do not repair yourself. Do not modify the device.
- Never use a damaged device. Replace defective or worn parts.
- Use the breast pump only for its intended use as described in these instructions for use.
- Do not use the breast pump while sleeping or overly drowsy.
- If tubing becomes moldy, discontinue use and replace tubing.

- Do not microwave or boil breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving (microwaving or boiling can also change the composition of breast milk).
- Clean and disinfect all parts that come into contact with your breast and breast milk prior to first use.
- Only use Medela recommended parts with your breast pump.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain. If you are unsure about correct breast shield sizing please visit **www.medela.com** or see a lactation consultant/breastfeeding specialist who can help you get a proper fit.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful. See **chapter 7.1** for more information.
- If pumping at high altitudes, including in an airplane, consider pumping more often or longer if you feel there is milk remaining in your breasts after your pumping session.

NOTICE


Take appropriate care in handling components:

- Plastic parts become brittle when frozen and may break when dropped.

- Parts may become damaged if mishandled (e.g., dropped, over-tightened, or knocked over).
- Do not use the breast milk from components that show signs of damage.

Contents

1. Important safety information	3	11. Maintenance and care	16
2. Intended purpose	6	11.1 Power adaptor care	16
2.1 Indications for use	6	11.2 Battery maintenance	16
2.2 Contraindications	6	11.3 Long-term storage	16
3. Product description	6	12. Troubleshooting	17
3.1 Description of the device	6	13. Ordering information	18
3.2 Your breast pump system features	6	14. Warranty	18
4. Getting started	7	15. Disposal	19
4.1 Cleaning before first use	7	16. Meaning of symbols	20
4.2 Charging the battery	7	17. International regulations	22
5. Disassembling your collection cups	7	17.1 Electromagnetic compatibility (EMC)	22
6. Cleaning	7	17.2 Radio transmission	24
6.1 Cleaning overview	8	18. Technical specifications	25
6.2 Washing	9		
6.3 Disinfecting	10		
6.4 Cleaning the motor unit	10		
7. Assembling your breast pump	11		
7.1 Selecting the correct breast shield size	11		
7.2 Assembling the collection cups	11		
7.3 Preparing for single pumping	12		
7.4 Preparing for double pumping	12		
7.5 Placing collection cups within your bra	12		
8. Operating your breast pump	12		
8.1 Mimicking the baby's rhythm (2-Phase Expression® technology)	12		
8.2 Buttons and status light	13		
8.3 Bluetooth® connectivity	13		
9. Expressing your breast milk	14		
10. Handling of breast milk	15		
10.1 Removing the collection cups from your bra	15		
10.2 Pouring milk from collection cups	15		
10.3 General storage guidelines	15		
10.4 Freezing	15		
10.5 Thawing	15		

EN Illustrated items or instruction steps are marked with the figure and detail number (e.g. )

2. Intended purpose

2.1 Indications for use

The Swing Maxi™ Hands-free breast pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts. The Swing Maxi™ Hands-free breast pump is intended for a single user. The breast pump is intended to be used in a home environment.

2.2 Contraindications

There are no known contraindications for the Swing Maxi™ Hands-free breast pump.

3. Product description

3.1 Description of the device

The Swing Maxi™ Hands-free breast pump is a personal-use electric breast pump that includes 2-Phase Expression® technology and is capable of single and double pumping. Mode of operation: Continuous. The operating life of the Swing Maxi™ Hands-free breast pump is 275 hours. Only the parts of the collection cups come into contact with the breast or milk. The motor unit generates an intermittent vacuum. For milk expression, the tubing conveys the resulting air flow via a membrane, which separates the milk from the pulsing air, to the breast.

Breast pump model number

101046709

The model number (REF) can be found near the UPC code on the package your breast pump came in.

Your breast pump system features

Refer to figure 1 on the foldout pages.

1-01 Hands-free Collection Cups

- 1-01-1 Breast shield 21mm (2pcs) (Applied part)
- 1-01-1 Breast shield 24mm (2pcs) (Applied part)
 - 1-01-2 Nipple Tunnel
- 1-01-3 O-ring (4pcs)
- 1-01-4 Yellow membrane (2pcs)
 - 1-01-5 Alignment Tab
 - 1-01-6 Valve
- 1-01-7 Outer shell (2pcs)
 - 1-01-8 Pouring spout
 - 1-01-9 Tubing connection port

1-02 Tubing

- 1-02-1 Short ends of tubing (connect to collection cups)
- 1-02-2 Long end of tubing (connects to motor unit)
- 1-02-3 Tubing holder (use for single pumping)

1-03 Motor unit (Applied part)

1-04 Power adaptor with cable

For ordering information about accessories see **chapter 13**.

All pictures shown in these instructions are for illustration purpose only. The actual product may vary due to local differences or product enhancement. Medela reserves the right to substitute any component or accessory with a replacement of equivalent performance.

4. Getting started

4.1 Cleaning before first use

It is important to do the following before using your breast pump for the first time:

1. **Disassemble and clean all parts** (except for the tubing) according to the instructions –see **chapter 5** and **6**.
2. **Please ensure to boil all parts of the collection cups before first use.** Do **not** clean or boil the tubing.
3. **When assembling, please ensure the cups are fully sealed until you hear a snap** – see **chapter 7.2**.

WARNING

Only use the power adaptor and cable that comes with the breast pump. An inadequate power source or cable may result in a fire hazard, electric shock or malfunction equipment due to overvoltage. Do not expose the breast pump or battery to excessive temperature.

4.2 Charging the battery

Charge the battery for **two hours**. Refer to figure  2 on the foldout pages.


1.  2-01 Connect the power cable to the motor unit.
 2.  2-02 Connect the power cable to the power adaptor.
 3.  2-03 Plug the power adaptor into a wall outlet.
→ During charging, the status light is pulsing in green color.
→ When charging is complete, it lights up steadily in green color.
- For details on the status light see **chapter 8.2**.

Information

- The warranty is only valid when using the Medela power adaptor to charge your pump (Model MS-V2000U050-010B0-* / MS-H2000U050-010B0-*).
- You can use the breast pump while charging the battery. However, it will charge in a slower mode when the pump is running.
*country code

5. Disassembling your collection cups

Refer to figure  3 on the foldout pages.

1.  3-01 Open the collection cup using the tabs (3-01-1) at the bottom.
2.  3-02 Remove the yellow membrane from the breast shield.

6. Cleaning

It is important you do the following before using your collection cup parts for the first time: Disassemble all parts and clean by following the cleaning instructions.

For additional cleaning guidelines see the Center for Disease Control website: <https://www.cdc.gov/healthywater/hygiene/healthychildcare/infantfeeding/breastpump.html> (U.S).

Supplies needed:

- Dish soap
- Clean soft brush
- Clean wash basin or dishwasher
- Drinking-quality water
- Clean pot for boiling water

Parts to wash or disinfect:

- Breast shields
- Outer shells
- Yellow membranes

CAUTION

Separate and wash all parts that are in contact with breast milk immediately after use. This will help to remove breast milk residue and prevent growth of bacteria.

- Wash hands thoroughly with soap and water.
- Only use drinking-quality tap water or bottled water for cleaning.
- Do not place collection cup parts directly in the sink to rinse or wash. Use a clean wash basin used only for infant items.
- Do not use a dish towel to rub or pat items dry.
- Do not use disinfectants for cleaning.

6.1 Cleaning overview

Washing and disinfection are two different activities. They must be done separately to protect you and your baby, and to maintain the performance of your breast pump.

Wash – To clean the surfaces of the parts by physically removing milk residues and other possible contaminants.

Disinfect – To efficiently lower the number of germs that may be present on the surfaces of the washed parts.

	Breast shield	Outer shell	Yellow membrane	O-ring
When to clean				
Before first use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Disinfect the parts according to chapter 6.3			During regular washing and disinfecting, you do not need to remove the O-ring from the breast shield.
After each use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2			
Once per day after use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Disinfect the parts according to chapter 6.3			
When to clean	Tubing		Motor unit	
As needed	Do not clean or rinse the tubing under any circumstance.		Wipe with clean, damp towel. Refer to chapter 6.4 .	

NOTICE

The Swing Maxi™ Hands-free breast pump has a closed system so milk cannot get into the tube when the pump is used as instructed. **Do not clean or rinse the tubing under any circumstance.**

• **Inspect the tubing before use. If you find condensation, breast milk or mold in the tubing, discontinue use and replace tubing.**

Refer to chapter 13 for ordering information.

• Do not run the breast pump with wet tubing; doing so will cause damage to your breast pump.


- Take care not to damage parts of the collection cups during cleaning.
- Store the dry collection cups in a clean bag/container until next use.
- Do not store wet or damp parts as mold may develop.
- If you notice a white residue on your collection cup parts after boiling, your water may have a high mineral content. Remove residue by wiping parts with a clean towel and allow to air dry.
- Distilled water is recommended when boiling parts to prevent substantial mineral build-up over time, which may compromise your parts.

6.2 Washing


Refer to figure 4 on the foldout pages.

Rinse and wash by hand (before first use and after each use)


Do not place the parts directly in the kitchen sink for rinsing and washing. Use a dedicated wash basin for infant items.

1. **Rinse** the disassembled parts, except for the tubing, with cold, clear drinking quality water (approx. 20 °C / 68 °F).
2. **Clean** these parts with plenty of warm, soapy water (approx. 30 °C / 86 °F). Use a commercially available dish soap and clean soft brush.
3. **Rinse** the parts with cold, clear drinking quality water for 10 to 15 seconds (approx. 20 °C / 68 °F).
4. **Drying:** Allow to air dry on a clean, unused dish towel or paper towel.

or

Wash in dishwasher (maximum once per day)

1. **Dishwashing:** Place the disassembled parts, except for the tubing, on the top rack or in the cutlery section. Use a commercially available dishwashing detergent.
2. **Rinse the yellow membrane** under running water after dishwashing.
3. **Drying:** Allow to air dry on a clean, unused dish towel or paper towel.


i Information

By using the dishwasher, parts may become discolored. This will not impact the functioning of the parts. **Never wash the tubing.** See **chapter 6.1** for more details.

6.3 Disinfecting

Refer to figure 5 on the foldout pages.

Disinfect on stovetop (before first use and once per day after use)

1. Cover the disassembled, washed parts, except for the tubing, with plenty of water and boil at least for 10 minutes. Ensure that the parts are not in direct contact with the bottom of the pot.
2. **Drying:** Allow to air dry on a clean, unused dish towel or paper towel. Do not store parts in an airtight container/bag if moist. It is important that all residual moisture dries.

or

Disinfect in microwave (once per day after use)

1. Use Quick Clean™ Micro-Steam™* bags in the microwave in accordance with the instructions on the bags. (*Refer to local website/shops for availability in your country.)
2. **Drying:** Allow to air dry on a clean, unused dish towel or paper towel. Do not store parts in an airtight container/bag if moist. It is important that all residual moisture dries.

10

6.4 Cleaning the motor unit

Supplies needed:

- Drinking-quality water
- Dish soap
- Clean dish towel or paper towel

Parts needed:

- Motor unit

⚠ CAUTION

To prevent electric shock and damage to the material:

- Unplug the breast pump before cleaning the motor unit.
- Do not immerse the motor unit in water; do not run water over the motor unit.

NOTICE

- Do not spray or pour liquid directly onto the pump.

1. Switch the breast pump off.
2. Unplug the breast pump from the power source.
3. Wipe the motor unit with a clean towel moistened with soapy water.
4. Wipe the motor unit dry with a clean, dry towel.

7. Assembling your breast pump

7.1 Selecting the correct breast shield size

Make sure to use the correct breast shield size. It is an essential component of effective pumping, helping to optimise milk flow.

Refer to figure  6 on the foldout pages.

1. Select the breast shield size by measuring the diameter of your nipple ( 6-01).
2. Select the breast shield size that is slightly larger than your nipple ( 6-02). 21mm and 24mm size shields are included in the pack.

For more information on breast shield sizing and other available sizes, visit www.medela.com or see **chapter 13** for ordering information.

7.2 Assembling the collection cups

Only use clean and dry parts for assembly. Refer to **chapter 6** for cleaning.

Parts needed:

- Tubing
- Breast shield with O-ring
- Motor unit
- Yellow membrane
- Outer shell

CAUTION

To prevent contamination of your milk:

- Wash hands thoroughly with soap and water before touching the motor unit, parts of the collection cups and breasts.
- Dry your hands with a fresh towel or a single-use paper towel.
- Avoid touching the inside of the collection cup components.

To make sure that your breast pump works properly and safely:

- Check collection cup components for wear or damage before use. Replace if necessary.
- Always inspect all parts prior to use for cleanliness. If dirty, clean according to **chapter 6**.
- Only use genuine Medela parts. See **chapter 13** for details.

NOTICE

To prevent damage to the breast pump all components must be completely dry before use.

Refer to figure  7 on the foldout pages.

1.  7-01 Wash hands thoroughly with soap and water before touching the motor unit, parts of the collection cups and breasts. Dry your hands with a clean towel.
2.  7-02 Check and confirm that the O-ring is applied to the breast shield.
3.  7-03 First assemble the valve of the yellow membrane on the front of the breast shield.
4.  7-04 Align the disc of the yellow membrane on the designated circular slot using the alignment tab at the top to ensure that it is at the right place. Press around the edges to make sure it is fully assembled.
5.  7-05 Insert the pouring spout into the opening on the shield.
5.  7-06 Push the outer shell and breast shield together tightly, ensuring it is fully sealed all the way around until you hear a snapping sound.
7.  7-07 Connect the long end of the tubing to the motor unit as far as it will go.
8.  7-08 Connect the short end of the tubing into the tubing connection port on the collection cup.


i Information

O-ring fits all breast shield sizes. To disassemble and reassemble the O-ring, follow the guidance below.

- **To disassemble the O-ring:** Start by pulling out the bottom tab (as shown in the left image below) from the breast shield.
- **To reassemble the O-ring:** Start from the top (as shown on the right image below) of the breast shield, stretch it over the outer edge of the breast shield. Ensure it goes on with the ridges facing outwards. Press along the groove, checking that it is in place all the way around.

To disassemble:


Pull out the bottom tab from here

To reassemble:

Start placing the O-ring from here.

Ensure that the ridges are straight

7.3 Preparing for single pumping

Refer to figure  on the foldout pages.

1. Assemble one collection cup as described in **chapter 7.2**.
 -  8-01 Insert the unused tubing end into the Y-shaped tubing holder.
 -  8-02 Correctly assembled system (for single pumping).
2. Place the collection cup in your bra according to **chapter 7.5**. Start pumping as described in **chapter 9**.

8. Operating your breast pump

8.1 Mimicking the baby's rhythm (2-Phase Expression® technology)

Babies use a “two-speed” sucking action. At first, they suck quickly, and then they switch to a slower, deeper sucking rhythm once milk starts to flow. The quicker sucking action stimulates the milk ejection reflex and triggers the milk flow; the slower sucking draws the milk out for the feeding phase.

The pump mimics the baby's natural sucking rhythm by using a “two-speed” pumping action, the so called “2-Phase Expression® technology”.

i Information

- Double pumping saves time and increases the nutritional value of the milk.

7.4 Preparing for double pumping

Refer to figure  on the foldout pages.

1. Assemble the two collection cups as described in **chapter 7.2**.
 -  9-01 Correctly assembled system (for double pumping).
2. Place the collection cups in your bra according to **chapter 7.5**. Start pumping as described in **chapter 9**.

7.5 Placing collection cups within your bra

Refer to figure  on the foldout pages.


1.  10-01 Ensure that the nipple is centered in the tunnel of the breast shield.
2.  10-02 The Hands-free Collection Cups can be worn under your regular nursing bra. Your bra will hold the collection cup in place. If you find it uncomfortable or painful during use, try switching to a looser, more stretchy nursing bra.
3.  10-03 The tubing can be connected facing inwards (a - towards the breasts) or outwards (b - towards the arms) based on personal preference. Ensure that the tubing is not pinched or bent during use.
 - Kinking could impact the vacuum performance.
4.  10-04 The Hands-free Collection Cups can be used in a vertical or slightly tilted position for better visibility of the nipple and milk ejection. The scale on the collection cup gives you an approximate indication of the volume during pumping.

The breast pump begins with a fast sucking action called the **stimulation phase**. When you see milk flow or feel milk ejection, you can then press the **“Let-down” button** to switch to the **expression phase**, which is slower than the previous one. If you do not press the Let-down button, the expression phase will automatically start after the pump has been running for one minute in the stimulation phase.

8.2 Buttons and status light

Refer to figure 11 on the foldout pages.

<ul style="list-style-type: none"> 11-01 Let-down button (switches between stimulation and expression phase) 11-02 Increase vacuum button 11-03 Decrease vacuum button 11-04 On/Pause/Off button 11-05 Status light 		<p>Status light: Mode of pump</p> <ul style="list-style-type: none"> Light is off Pump is off Light is flashing Pump is on, stimulation phase Lights up steadily Pump is on, expression phase Light is pulsing Pump is on, paused 	<p>Status light: Battery level</p> <ul style="list-style-type: none">  Good  Low, approximately 10 minutes left  Very low, charge the battery immediately  Light is pulsing: Charging Lights up steadily: Battery fully charged
--	---	--	--

8.3 Bluetooth® connectivity

Bluetooth® is a wireless technology for exchanging data over short distances. This is used to connect your Swing Maxi™ Hands-free breast pump to the Medela Family™ smartphone app. By connecting your breast pump to the Medela Family™ app you can automatically transfer your pumping data (session length, phases & levels) and manually input the amount of milk expressed. To benefit from the features of the app you need to create an account in the Medela Family™ app.

Installing the Medela Family™ app

1. Download and install the Medela Family™ app to your mobile device. The app is available free of charge for Apple® iOS and Android™ operating systems.
2. When starting the Medela Family™ app for the first time you are asked to give the app permission to connect to Bluetooth® devices. Please give the app permission to allow the app to connect to your Medela pump.
 - >Android™ phones only: You will also be asked to give the app permission to use your location. We don't track the position of your phone, but this is a requirement of the Android™ operating system to connect to Bluetooth® devices.
 - >If you denied the permissions and now you can't connect to the pump:
 1. Open the Settings app of your smartphone, then (depending on your smartphone brand, you might need to select a submenu called "Apps") scroll down and find Medela Family™.
 2. There, make sure Bluetooth® permission is authorised. On Android™: also authorise location permission.
3. Open the Medela Family™ app on your mobile device and register by means of the app.

Connecting (pairing) your Swing Maxi™ Hands-free breast pump

1. Make sure that Bluetooth® is enabled on your mobile device.
2. In Medela Family™, click on "More", navigate to "Pump settings" and follow the set-up and pairing instructions on the screen, which include the following steps:
3. Press and hold the Let-down button  for two seconds to start the pairing procedure.
 - The status light is pulsing in blue color to show that the breast pump is ready for pairing with your Bluetooth® device.
4. Successful pairing will be indicated with the blue light being on steadily for two seconds.
5. After initial setup, the status light will turn off.
6. When reconnecting to Medela Family™, the status light changes to blue color for two seconds.

Saving your pumping session into the app

After you finished a session press the On/Pause/Off button of the breast pump for 2 seconds to transmit the data to the Medela Family™ app and to switch off the breast pump.

Removing your Swing Maxi™ Hands-free breast pump from your smartphone

1. Inside the Medela Family™ app go to "More" and then to "Pump settings".
2. Tap on the pump you want to remove and then tap on "Remove this pump" button.

i Information

- The Swing Maxi™ Hands-free breast pump has the ability to store your last 30 pump sessions in its memory. If the pump is not connected prior to session 31 your history will be overwritten with new session data.
- When connected to the pump, your Medela Family™ account will log an infinite number of historical sessions.
- Swing Maxi™ Hands-free uses Bluetooth® technology. Please follow airline rules for the use of portable electronic devices when flying.
- Removing the Swing Maxi™ Hands-free breast pump from your mobile device will not delete pumping information stored in the Medela Family™ app.
- If you reconnect your breast pump and mobile device the last 30 records stored within your pump will be transferred.

To learn more, please visit www.medela.com/app.

Trademark recognition: The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Medela AG is under license. Apple is a trademark of Apple Inc., registered in the U.S. and other countries. Android is a trademark of Google LLC.

9. Expressing your breast milk


⚠ CAUTION

- Do not try to express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful.

NOTICE

Make sure that the tubing is not kinked or pinched while pumping.

Refer to figure 12 on the foldout pages.


1. Make sure that the collection cups are connected and positioned properly.
2. 12-01 Press the On/Pause/Off button  to start pumping.
→ The stimulation phase begins and the status light is flashing.
3. 12-02 Adjust the vacuum with the Increase vacuum  and Decrease vacuum  buttons to find a comfortable level.
4. 12-03 As soon as your milk begins to flow, press the Let-down button 
→ The expression phase begins and the status light is on steadily.
→  If the Let-down button is not pressed within one minute, the breast pump will automatically go to expression phase.
5. 12-04 Adjust the vacuum with the Increase vacuum  and Decrease vacuum  buttons to find your maximum comfort vacuum (see information below).

Once you are pumping in the **expression phase**, increase the vacuum level until pumping feels slightly uncomfortable (not painful), then decrease the vacuum by 1 level with the Decrease vacuum  button.

6. When your pumping session is over, press the On/Pause/Off button  to stop the pumping operation.
→ The breast pump stops automatically after 30 minutes if there is no user interaction with the pump.
7. To switch the breast pump off press and hold the On/Pause/Off button  for two seconds.

i Information

- **Maximum comfort vacuum** is the highest vacuum setting where pumping still feels comfortable. This is different for every mother.
- Stimulation should be at a comfortable vacuum level, striving for a maximum level is not necessary.
- Reassess your maximum comfort vacuum throughout your pumping experience. It may change throughout each stage of lactation.
- The pump is not designed to be used while lying down or sleeping. However, you can slightly lean forward or backward during your pumping session.
- The maximum volume/capacity for the collection cup is 150ml/5oz. When this limit is exceeded, transfer the milk to a bottle before continuing to express.

Pausing

In both stimulation and expression phases the breast pump can be paused. This allows you to readjust yourself or your collection cups, or attend to your surroundings. For pausing, press the On/Pause/Off button  while the breast pump is running.

→ The pumping operation stops and the status light is pulsing.
To resume pumping, press the On/Pause/Off button  again.

i Information

- If pumping is not resumed within two minutes, the breast pump switches off automatically.
- To resume pumping, you will need to start a new session.
- If you need to pause for more than two minutes, we recommend restarting with stimulation phase.
- Do not forget to record milk output in the Medela Family™ app.

10. Handling of breast milk

10.1 Removing the collection cups from your bra

Refer to figure  13 on the foldout pages.

1.  13-01 While the pump is still running, slightly lean forward and ensure all milk from the nipple tunnel enters the collection container. Turn off the pump.
2.  13-02 Detach the tubing from the collection cups.
3.  13-03 Carefully remove the collection cups from your bra.

10.2 Pouring milk from collection cups

1.  13-04 Place the collection cups in an **upright position** on a surface where they will not be easily knocked over. You can use the volume scale markings on the shield to read the collected milk volume. The scale on the outer shell gives you an approximate indication of the volume during pumping.
2.  13-05 **With the rounded outer shell facing down**, pour breast milk through the pouring spout into a bottle or container for storage. Do not use collection cups to store milk. Do not store assembled cups in the refrigerator.
3. Fully disassemble and clean after each use and disinfect once per day as described in **chapter 6**.

10.3 General storage guidelines

Storage guidelines for freshly expressed breast milk (for healthy term babies)

Room temperature 16 to 25°C (60 to 77°F)	Refrigerator 4°C (39°F) or colder	Freezer -18°C (0°F) or colder	Breast milk thawed in the refrigerator
Up to 4 hours is best	Up to 3 days is best	Up to 6 months is best	At room temperature: Up to 2 hours
*Up to 6 hours for milk expressed under very clean conditions	*Up to 5 days for milk expressed under very clean conditions	*Up to 9 months for milk expressed under very clean conditions	Refrigerator: Up to 24 hours Do not refreeze!

*Guidelines for expressing milk under very clean conditions: Before expressing breast milk, mothers should wash their hands with soap and water or an alcohol-based hand sanitizer. The pump parts and pumping area must be clean. Breasts and nipples do not need to be washed before pumping.

- These guidelines for storage and thawing of breast milk are a recommendation. Contact your lactation consultant or breastfeeding specialist for further information.
- Store breast milk in the coldest spot of the refrigerator (at the back of the glass shelf above the vegetable compartment).
- For information on storing and thawing your breast milk, visit www.medela.com or visit www.breastmilkguidelines.com (U.S.) or <https://www.womenshealth.gov/breastfeeding/pumping-and-storing-breastmilk> (U.S.) or <https://www.medela.ca/breastfeeding/moms-journey/storing-and-thawing-breast-milk> (Canada).

10.4 Freezing

NOTICE

- Repeated freeze/thaw cycles destroy the structure of the milk. Therefore, never refreeze breast milk.
- Freeze expressed breast milk in food-grade bottles or milk storage bags. Do not fill the bottles or bags more than 3/4 full to allow space for possible expansion.
- Label the bottles or bags with the date and volume of expression.
- Plastic bottles and parts become brittle when frozen and may break when dropped.
- Do not use the breast milk from bottles or components that show signs of damage.

10.5 Thawing

CAUTION

Do not thaw or warm breast milk in a microwave or a pan of boiling water. This helps preserve important components and prevents burns.

- To preserve breast milk components, thaw the milk overnight in the refrigerator.
- Alternatively, hold the bottle or bag under warm water (max. 37°C / 98.6°F) and use immediately, as soon as it has been thawed.

NOTICE


Gently swirl the bottle or bag to mix any fat that has separated. Avoid shaking or stirring the milk.

11. Maintenance and care

11.1 Power adaptor care

WARNING

If the power adaptor outer casings become loose or are separating, stop use and contact Medela immediately. Separation of the casing can result in electric shock.

NOTICE

To keep your power adaptor and cable in good working condition, pay attention to the following: Do not unplug the cable or the power adaptor by pulling on the cord, but hold it by the connector housing and adaptor body.

Traveling internationally

Information

When traveling internationally, we recommend purchasing a travel adaptor (not included) to match the corresponding wall outlet. Please consult with the country you are visiting to find out what type of adaptor will work best. The power adaptor of the breast pump is two-pronged and does not require a grounded outlet. This breast pump has a lithium-ion battery.

Information

Carrying such batteries on airplanes may be restricted by the country you are visiting. Please consult with the country you are visiting to find out if there are any restrictions that pertain to traveling with lithium-ion batteries.

11.2 Battery maintenance

Your breast pump contains a built-in (not replaceable) lithium-ion rechargeable battery.

NOTICE

- In case of a completely discharged battery the breast pump cannot be used instantly after connecting it to the power adaptor. Allow the battery to be charged for at least 15 minutes before attempting to switch the pump on.
 - For optimum pump and battery performance use the power adaptor that comes with the breast pump (Model MS-V2000U050-010B0-* / MS-H2000U050-010B0-*).
 - Make sure that the voltage of the power adaptor is compatible with the power source.
 - Be aware of the risk of USB-C-socket of the pump or battery being shortcircuited by metal objects or fluids.
- *country code

Information


To preserve the battery life cycle:

- Keep your breast pump in a cool place.
- Recharge the battery before it completely drains. This is better for the life cycle of the battery
- Fully charge the battery before storing your breast pump. Refer to **chapter 8.2** for information on the battery charging state.
- If your breast pump has been stored in a hot location, it may not run on battery power right away. To resume normal battery function, allow the breast pump to cool down for one hour. During this time, you can power your breast pump with the power adaptor plugged into a wall outlet.

11.3 Long-term storage

Your breast pump contains a rechargeable battery. To preserve your breast pump, when not in use for an extended amount of time, it will automatically transition into storage mode. To remove from storage mode, plug your breast pump into an external power source and charge the battery for two hours. Do not store collection cups in assembled state over an extended period of time.


Before storing your breast pump for an extended amount of time remember to:

Clean the collection cup (refer to **chapter 6**)

Clean the motor unit (refer to **chapter 6.4**)

Make sure that the battery is charged (refer to the following information)

It is best to fully charge the battery before storing your breast pump for an extended period of time. Refer to **chapter 8.2** for information on the battery charging state.

12. Troubleshooting

In case of an unexpected behavior of your breast pump check with the troubleshooting table if you find the issue in the column “Problem” and follow the instructions in the column “Solution”.

Problem	Solution	Problem	Solution
The breast pump is not working after you pressed the On/Pause/Off button	<ul style="list-style-type: none"> • Charge for 15 minutes before restarting the pump. Charging is indicated by the status light pulsing in green. • If the pump is still not working contact Medela Customer Service. 	There is low or no suction	<ul style="list-style-type: none"> • While pumping, make sure the breast shields form a complete seal around the breast. • When single pumping, make sure that the unused tubing end is correctly plugged into the tubing holder. • If suction does not improve after having followed these steps, contact Medela Customer Service.
The breast pump generates only constant vacuum instead of sucking cycles	<ul style="list-style-type: none"> • Contact Medela Customer Service. 	The breast pump is not responding as expected	<ul style="list-style-type: none"> • Press and hold the On/Pause/Off button for two seconds to turn the breast pump off, then restart the breast pump. • If there is no change, contact Medela Customer Service.
The status light is flashing rapidly in red when you try to pump or charge the breast pump	<ul style="list-style-type: none"> • Charging/pumping is not possible when the battery temperature is too high. Allow the motor unit to cool down, then resume charging/pumping. • If the status light is still flashing rapidly in red, a device fault occurred. Contact Medela Customer Service. 	The motor unit's exterior got wet	<ul style="list-style-type: none"> • Unplug the motor unit from the power source and turn it off. • Dry off the exterior of the motor unit.
There is low or no suction	<ul style="list-style-type: none"> • Make sure that all pump set components are clean and dry and that connections are secure. 		

Problem	Solution
The motor unit has been submerged in water	<ul style="list-style-type: none"> • Unplug the motor unit from the power source and turn it off. • Contact Medela Customer Service.
Data cannot be sent	<ul style="list-style-type: none"> • The Bluetooth® function “bluetooth picture” on the mobile device is turned off. Turn the Bluetooth® function on and try to send the data again. • Pairing between the two devices has not been completed. Perform pairing instructions again. Refer to chapter 8.3. • The mobile device is too far away from the breast pump. Move the device and breast pump closer together. • The application on the mobile device is not ready. Check the application.

If you have not resolved the problem with your breast pump or you have further questions, please contact Medela Customer Service. For contact data visit www.medela.com/contact-us. Under “Country” choose your country.

14. Warranty

This product is warranted by Medela to the original retail purchaser to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and detachable components) from the date of purchase. Warranty can only be claimed in the country of purchase. In the event of a defect, Medela will repair or, at Medela’s option, replace this product, without charge for such replacement, parts or labor. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

13. Ordering information

Accessories

Article number	Product
Depending on your region, check local website/shops for availability.	21mm Breast shields compatible with Medela Hands-free Collection Cups*
	24mm Breast shields compatible with Medela Hands-free Collection Cups*
	27mm Breast shields compatible with Medela Hands-free Collection Cups*
	O-ring compatible with Medela Hands-free Collection Cups
	Membrane compatible with Medela Hands-free Collection Cups
	Y-shaped tubing compatible with Swing Maxi™ Hands-free
	Standalone Set Hands-free Collection Cups

* If the result of expression is unsatisfactory or expression is painful, contact your lactation consultant or breastfeeding specialist. A different size of breast shield may help make expression more comfortable and effective. In case of difficulties in finding your desired accessory, please contact Medela Customer Service.

For contact data visit www.medela.com/contact-us. Under “Country” choose your country. For replacement of lost or defective parts please contact Medela Customer Service. Other Medela products are available at www.medela.com.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO A DURATION OF 1 YEAR FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE OR IN YOUR COUNTRY.

Before you make a claim under this warranty, it may save time and expense to call Medela Customer Service (see contacts below). You may also call this number for additional information concerning this warranty. If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, prepaid, together with your dated bill of sale or other proof of purchase and a brief statement of the problem to the following address:

USA
Medela LLC – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS
Medela Customer Service (toll free):
1-800-435-8316

Medela Canada
4160 Sladeview Crescent Unit #8
Mississauga, ON L5L 0A1
ATTENTION: RETURNS
info@medela.ca

Contact Customer Service for a return authorization number.
Returns are not accepted without an authorization number.

15. Disposal

Disposal of your breast pump

At the end of its operating life, separate the parts of your breast pump and dispose of according to the following instructions.

This appliance contains batteries that are non-replaceable.

When the battery is at end of life, the appliance shall be properly disposed of.

The integrated battery must not be removed for disposal.

If the battery should leak, the breast pump must be disposed of.

Breast shields, outer shells, tubing and membranes

The parts are made of plastics that are not harmful to the environment when disposed of as household waste. Recycle or dispose of according to local regulations.

Motor unit and power adaptor

Do not dispose of electric or electronic equipment together with unsorted municipal waste, but collect it separately.

In the European Union the manufacturer or its vendor must take back waste equipment. Other countries may have similar collection and recycling systems. Dispose of the waste equipment according to local regulations.

Inquire at the point of sale or contact your local authority for appropriate collection points for waste equipment.



The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment.

CAUTION

Hazardous materials present in electric/electronic equipment may harm human health and the environment if the waste equipment is not disposed of properly.

16. Meaning of symbols

The following tables explain the meaning of the symbols found on the product parts and its packaging.

	General safety alert symbol, points to information related to safety. ¹		Defines a temperature range (e.g. for operation, transport or storage). ²		Defines a relative humidity range (e.g. for operation, transport or storage). ³
	Defines an atmospheric pressure range (e.g. for operation, transport or storage). ⁴		Identifies the manufacturer. ⁵		Indicates the date of manufacturing. ⁶
	Indicates the serial number of the device. ⁷	IP20 IP22	Indicates the degree of protection against ingress of foreign objects and moisture. ⁸		Indicates compliance with international requirements for protection from electric shock (Type BF applied parts). ⁹
	Do not dispose of electric/electronic devices together with unsorted municipal waste (dispose of the device in accordance with local regulations). ¹⁰		Read and follow the instructions for use. ¹¹		The CE mark indicates conformity with the applicable European directives. ^{*, †}
	Indicates direct current. ¹²		Indicates the location of the On/Pause/Off button. ¹³		Indicates alternating current. ¹⁴
	Indicates the polarity of the USB power output. ¹⁵		Indicates that the device is a Class II electrical appliance (double insulated). ¹⁶		Indicates that the device is for indoor use only. ¹⁷
	Indicates that a specific separate power adaptor is required for connecting the device to the supply mains and specifies its model/type MS-V2000U050-010B0-*/MS-H2000U050-010B0-**. ¹⁸		Indicates that the appliance is intended to be usable up to the maximum altitude of 5000 m. [*]		Indicates the efficiency level for external power supplies. ¹⁸
	Indicates the part number of the product. ¹⁹		Contains fragile goods. Handle with care. ²⁰		Keep away from sunlight. ²¹
	Keep away from rain. Keep in dry conditions. ²²		Indicates that an item can technically be recycled. ²³		Indicates that the device is a short-circuit proof safety isolating transformer. ²⁴
	Indicates the compliance with the requirements of the Federal Communications Commission. ^{25, †}		Indicates the material from which an item is made. ^{26, †}		Indicates the model number of the device. ²⁷
	Indicates that this device contains Bluetooth® wireless technology. [†]		Indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment. ^{28, †}		Indicates that the electronic device incorporating transformer(s) and electronic circuitry(ies), converts electrical power into single or multiple power outputs. ²⁹

The UL LISTED mark indicates that the product is manufactured in compliance with UL safety requirements for USA and Canada. †

The NOM (Norma Oficial Mexicana) NYCE mark indicates that the product is found to be in compliance with the applicable Mexican requirements.†

Indicates the lot/batch number.³⁰

Indicates the compliance with PSE requirements (Japan, Product Safety Electrical Appliance and Material Safety Law).*

The BSMI (Bureau of Standard, Metrology and Inspection) logo indicates that the product is certified to meet Taiwan's requirements for safety and electromagnetic compatibility.*

Note: Symbols may vary depending on the market the unit is sold.

References

- 1 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.2 Symbol 2 General warning sign
- 2 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.7 Temperature Limit / ISO 7000-0632, Graphical symbols for use on equipment, Temperature Limit
- 3 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.8 Humidity Limit / ISO 7000-2620, Graphical symbols for use on equipment, Humidity Limit
- 4 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.9 Pressure Limit / ISO 7000-2621, Graphical symbols for use on equipment, Pressure Limit
- 5 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.1.1 Manufacturer
- 6 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.1.3 Manufacturing Date / ISO 7000-2497, Graphical symbols for use on equipment, Date of manufacture
- 7 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.1.7 Serial number / ISO 7000-2498, Graphical symbols for use on equipment, Serial number
- 8 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.3 Symbol 2 IP Code
- 9 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 20 Type BF applied parts
- 10 EN 50419, Marking of Electrical and Electronic Equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE).
- 11 IEC 60601-1, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance, Table D.2 Symbol 10 Refer to instruction manual/booklet
- 12 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 4 Direct current
- 13 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 29 Stand-by
- 14 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 1 Alternating current
- 15 IEC 60950-1 Information technology equipment - Safety - Part 1: General requirements, Clause 1.7.1 / IEC 60417 Graphical symbols for use on equipment
- 16 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 9 Class II equipment
- 17 IEC 60417-5957, Graphical symbols for use on equipment, For Indoor use only
- 18 US Department of Energy Requirement, 10 CFR Part 430
- 19 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.1.6 Catalogue number / ISO 7000- 2493, Graphical symbols for use on equipment, Catalogue number
- 20 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.1, Fragile, handle with care / ISO 7000-0621, Graphical symbols for use on equipment, Fragile, handle with care
- 21 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.2 Keep away from sunlight / ISO 7000-0624, Graphical symbols for use on equipment, Keep away from sunlight
- 22 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.3.4, Keep dry / ISO 7000-0626, Graphical symbols for use on equipment, Keep dry
- 23 ISO 7000-1135, Graphical symbols for use on equipment, General symbol for recovery/recyclable
- 24 IEC 60417-5947, Graphical Symbols for Use on Equipment, Safety isolating transformer, short-circuit proof
- 25 Code of Federal Regulations, Title 47, Part 15b / 15c.
- 26 97/129/EC, Commission Decision establishing the identification system for packaging materials.
- 27 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.1.10 Model number / IEC 60417-6050 Graphical Symbols for Use on Equipment, Model number
- 28 TÜV (Technical Inspection Association) mark indicates that the product is manufactured in compliance with UL safety requirements for USA and Canada (USA: UL60950-1, CAN: CSA C22.2 NO. 60950-1).
- 29 IEC 60417-6190, Graphical Symbols for Use on Equipment, Power supply unit, switch mode; switch mode power supply unit; SMPS
- 30 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.5 Batch code / ISO 7000-2492, Graphical symbols for use on equipment, Batch code.

* These symbols are not applicable for the US, Canada, Mexico market.

† These symbols are not derived from standards.

17. International regulations

17.1 Electromagnetic compatibility (EMC)

The Swing Maxi™ Hands-free breast pump is intended to be used in a home environment and is EMC-tested in conformity with the requirements of the standard IEC 60601-1-2.

The breast pump needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in these instructions for use. Portable and mobile RF communications can affect the breast pump.

WARNING

To prevent malfunction resulting from electromagnetic interference:

- The electric breast pump Swing Maxi should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the electric breast pump Swing Maxi should be observed to verify normal operation in the configuration in which it will be used.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies, RFID can affect the electric breast pump and should be kept at a distance of at least 30 cm away from the device.
- Use of accessories or cables other than those provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

NOTICE

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

NOTICE

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Guidance and manufacturer's declaration – electromagnetic emissions

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Emission tests	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	
Harmonic emissions IEC 61000-3-2	Class A (only valid with power adaptor supplied by Medela)	The breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliant	

Guidance and manufacturer's declaration – electromagnetic immunity

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment. The Swing Maxi™ Hands-free breast pump has no essential performance but was tested for immunity to electromagnetic disturbances and passed using the following criteria:

1. No change in mode and operation of the breast pump.
2. The pump can be controlled by changing the vacuum level and/or pressing the On/Pause/Off button.

Reset (The following sequence can be used to reset the device in the event of an EM Disturbance)

1. Unplug the power adaptor from the wall outlet.
2. Unplug the power cable from the pump unit.
3. Attach a pump set to the tubing port (if none is attached).
4. Plug the power cable into the supply socket of the pump unit.
5. Plug the power adaptor back into the wall outlet.
6. Press the On/Pause/Off button to begin pumping.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact discharge +/- 15 kV air discharge	+/- 8 kV contact discharge +/- 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 5%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV, ± 1 kV Line-to-line	± 0,5 kV, ± 1 kV Line-to-line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 60601-1-2 IEC 61000-4-11	0 % U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U _T for 1 cycle at 0° 70 % U _T for 25 cycles at 50 Hz at 0° 0 % U _T for 250 cycles at 50 Hz at any one phase	0 % U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U _T for 1 cycle at 0° 70 % U _T for 25 cycles at 50 Hz at 0° 0 % U _T for 250 cycles at 50 Hz at any one phase	Mains power quality should be that of a typical commercial or hospital environment. If the user of the breast pump requires continued operation during power mains interruptions, it is recommended that the breast pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	It may be necessary to position the breast pump further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table of frequencies of portable and mobile transmitters for which the recommended separation distance is 30 cm (12 inches):

Band (MHz)	Service
380 - 390	TETRA 400
430 - 470	GMRS 460, FRS 460
704 - 787	LTE Band 13, 17
800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5
1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS
2 400 - 2 570	Bluetooth®, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7
5 100 - 5 800	WLAN 802.11 a/n

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the breast pump including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

17.2 Radio transmission

FCC (USA) statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

ISED (Canada) statement

This device complies with Industry Canada's licence-exempt RSSs. Operation is subject to the following two conditions:

- (1) this device may not cause interference; and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

18. Technical specifications

Summary of important technical specifications

Vacuum performance (at 500 m above sea level [954.62 hPa ambient pressure])
-45 to -245 mmHg (-60 to -327 hPa)
45 to 111 cpm

Power adaptor (Model: MS-V2000U050-010B0-* / MS-H2000U050-010B0-*)
*country code

Input	100 to 240 V~ 50/60 Hz max. 0.5 A
Output	5 V $\overline{\text{---}}$ 2 A 10 W
Average efficiency	78.7 %
No-load power consumption	less than 0.075W

Motor unit


Battery capacity and type: 3.6 V, 2500 mAh (nominal) Li-ion.
Size: 126 x 125 x 58 mm. Weight: 325 g. Ingress protection level: IP22.

Transmission standard: Bluetooth® version 4.2 BLE
Frequency range: 2402-2480 MHz
Max. power output: +4 dBm (class 2)

Materials touching skin or coming in contact with milk

- Breast shield : Polypropylene
- Outer shell : Polypropylene
- Membrane : Silicone
- O-ring : Silicone
- Tubing : Polyvinyl Chloride, Polycarbonate, Polybutylene Terephthalate

Not intentionally made of BPA (Bisphenol A).

 +35°C
+5°C
Operating temperature range
(5 °C to 35 °C)
(41 °F to 95 °F)

 93%
15%
Operating humidity range
(15% to 93% relative humidity)

 106kPa
70kPa
Operating atmospheric pressure range:
70kPa...106kPa

 +50°C
-20°C
Transport/storage
temperature range
(-20°C to 50°C)
(-4 °F to 122 °F)

 93%
15%
Transport/storage
humidity range
(15% to 93% relative humidity)

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