Hands-free Collection Cups

Collecteurs
Mains-libres

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Instructions
for use

EN ENGLISH
4-16

Mode
d'emploi

FR FRANÇAIS
17-27

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ES ESPAÑOL
28-38

Contents

1. Quick start	4
2. Important safety information	7
3. Intended use	8
4. Product description	8
5. Disassembling your Hands-free Collection Cups	8
6. Cleaning	8
7. Assembling your Hands-free Collection Cups	11
8. Handling of breast milk	13
9. Troubleshooting	14
10. Ordering information	14
11. Warranty	15
12. Disposal	15
13. Meaning of symbols	15
14. Technical specifications	16

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1-02

1-02-1 1-02-2

- 1-01 Hands-free Collection Cups
 - 1-01-1 Breast shield 21 mm (2pcs) (Applied part)
 - 1-01-1 Breast shield 24 mm (2pcs) (Applied part)
 - 1-01-2 Nipple tunnel
 - 1-01-3 O-ring (2pcs)
 - 1-01-4 Yellow membrane (2pcs)
 - 1-01-5 Alignment tab
 - 1-01-6 Valve
 - 1-01-7 Outer shell (2pcs)
 - 1-01-8 Pouring spout
 - 1-01-9 Tubing connection part
- 1-02 Tubing (2pcs)
 - 1-02-1 The flat ends of tubing (connect to collection cups)
 - 1-02-2 The round ends of tubing (connect to your personal use pump's tubing)

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







2. Important safety information

 **Before using this product for pumping, please read your Medela compatible personal use pump Instructions for Use for additional important safety information.**

Warnings identify all Instructions for Use that are important for safety. Failure to observe these Instructions for Use can lead to personal injury or damage to the product. The following symbols and signal words show the significance of the warnings:

WARNING

Can lead to serious injury or death.

CAUTION

Can lead to minor injury.

NOTICE

Can lead to material damage (not related to personal injury).

Information

Useful or important information that is not related to safety.

WARNING

To avoid health risks and reduce the risk of severe injury:

- The Hands-free Collection Cups are a single-user product. Use by more than one person may present a health risk and voids the warranty.
- Wash all parts that come into contact with your breast and breast milk after every use.
- The tubing poses a risk of strangulation when played with and small parts may be swallowed by toddlers. Close supervision is necessary when the Hands-free Collection Cups are used in the vicinity of children.

CAUTION

To avoid fire or burns:

- The Hands-free Collection Cups and accessories are not heat-resistant. Keep away from heated surfaces or open flames.
- Do not use near flammable materials.

To avoid health risks and reduce the risk of injury:

- Before each use visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- This product cannot be serviced or repaired. Do not repair yourself. Do not modify the device.
- Never use a damaged device. Replace defective or worn parts.
- Use the Hands-free Collection Cups only for their intended use as described in these Instructions for Use.
- If tubing becomes moldy, discontinue use and replace tubing.
- Do not microwave or boil breast milk. Microwaving can cause severe burns to a baby's mouth from hot spots that develop in the milk during microwaving. Microwaving or boiling can also change the composition of breast milk.

CAUTION

- Clean and sanitize all parts that come into contact with your breast and breast milk prior to first use.
- Use only with Medela compatible personal use pumps.
- If pumping is uncomfortable or causes pain, turn the pump unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- Choose the correct breast shield size for pumping. Using an unsuitable breast shield size may cause pain, discomfort and reduction in milk flow. Visit www.medela.com/fittingguide or see a lactation consultant/breastfeeding specialist who can assist with fit or if you can express only minimal or no milk or if expression is painful. See **chapter 7.1** for more information.

3. Intended use

3.1.1 Indications for use

The Hands-free Collection Cups reusable set is intended to be used by lactating women for expressing breast milk with the Medela double electric personal use breast pumps (Freestyle™, Pump in Style®, Swing Maxi™ / Duo). The pump set is a single-user product.

3.1.2 Contraindications

There are no known contraindications for the Hands-free Collection Cups.

4. Product description

Refer to figure  1 on Chapter 1.

For ordering information about accessories see **chapter 10**.

All pictures shown in these Instructions for Use are for illustration purpose only. The actual product may vary due to local differences or product enhancement. Medela reserves the right to substitute any component or accessory with a replacement of equivalent performance.

5. Disassembling your Hands-free Collection Cups

Refer to figure  2 on Chapter 1.

1.  2-01 Open the collection cup using the tab at the bottom.
2.  2-02 Remove the yellow membrane from the breast shield.

6. Cleaning

It is important you do the following before using your Hands-free Collection Cups for the first time: Disassemble all parts & clean following cleaning instructions.

For additional cleaning guidelines see the Center for Disease Control website:

<https://www.cdc.gov/healthywater/hygiene/healthychildcare/infantfeeding/breastpump.html> (U.S).

Supplies needed:

- Dish washing liquid
- Clean soft brush
- Clean wash basin or dishwasher
- Drinking-quality water
- Clean pot for boiling water

Parts to wash and sanitize:

- Breast shields (with O-ring)
- Outer shells
- Yellow membranes

CAUTION

Separate and wash all parts that are exposed to breast milk immediately after use. This will help to remove breast milk residue and prevent growth of bacteria.


- Wash hands thoroughly with soap and water.
- Only use drinking-quality tap water or bottled water for cleaning.
- Do not place the Hands-free Collection Cups' parts directly in a sink to rinse or wash. Use a clean wash basin used only for infant items.
- Do not use a dish towel to rub or pat items dry.
- Do not use disinfectants for cleaning.

NOTICE

Medela personal use breast pump has a closed system so milk cannot get into the tubing when the pump is used as instructed.

Do not clean or rinse the tubing under any circumstance.

- **Inspect the tubing before use. If you find condensation, breast milk or mold in the tubing, discontinue use and replace tubing. Refer to chapter 10 for ordering information.**
- Do not run your breast pump with wet tubing; doing so will cause damage to your breast pump.

NOTICE

Pay attention to the following:

- Take care not to damage parts of the Hands-free Collection Cups during cleaning.
- Store the dry Hands-free Collection Cups in a clean bag/container until next use.
- Do not store wet or damp parts as mold may develop.
- If you notice a white residue on your parts after boiling, your water may have a high mineral content. Remove residue by wiping parts with a clean towel and allow to air dry.
- Distilled water is recommended when boiling parts to prevent substantial mineral build-up over time, which may compromise your parts.

6.1 Cleaning overview

Washing and sanitizing are two different activities. They must be done separately to protect you and your baby, and to maintain the performance of your breast pump.

Wash – To clean the surfaces of the parts by physically removing milk residues and other possible contaminants.

Sanitize – To efficiently lower the number of germs that may be present on the surfaces of the washed parts.

Do not clean or rinse the tubing under any circumstance.

Tubing

When to clean	Breast shield 	Outer shell 	Yellow membrane 	O-ring 
Before first use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Sanitize the parts according to chapter 6.3			During regular washing and sanitizing, you do not need to remove the O-ring from the breast shield.
After each use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2			
Once per day after use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Sanitize the parts according to chapter 6.3			

6.2 Washing

Refer to figure  3 on Chapter 1.

Rinse and wash by hand (before first use and after each use)

Do not place the parts directly in a kitchen sink for rinsing and washing. Use a dedicated wash basin for infant items.

1. Rinse the disassembled parts, except for the tubing, with cold clear drinking-quality water (approx. 68 °F/20 °C).
2. Clean the parts with plenty of warm, soapy water (approx. 86 °F/30 °C). Use a commercially available dish soap and a clean soft brush.
3. Rinse the parts with cold clear drinking-quality water for 10 to 15 seconds (approx. 68 °F/20 °C).
4. Allow to completely dry after washing.

or

Wash in dishwasher (maximum once per day)

1. Place the disassembled parts, except for the tubing, on the top rack or in the cutlery section. Use a commercially available dishwashing detergent.
2. Rinse the yellow membrane under running water after dishwashing.
3. Allow to completely dry after washing.

Information

If using a dishwasher, parts may become discolored. This will not impact the functioning of the parts.

Drying

1. Allow to air-dry on a clean, unused dish towel or single-use paper towel.

6.3 Sanitizing

Refer to figure 4 on Chapter 1.

Sanitize on stovetop (before first use and once per day after use)

1. Cover the disassembled, washed parts, with plenty of water and boil at least for 10 minutes. Ensure that the parts are not in direct contact with the bottom of the pot.
2. Allow to dry after sanitizing.

or

Sanitize in microwave (once per day after use)

1. Use the Medela's Quick Clean Microwave bags* bags in a microwave in accordance with the Instructions for Use on the bags.

2. Allow to dry after sanitizing.

* Refer to local website/shops for availability in your country.

Drying

1. Allow to air-dry on a clean, unused dish towel or single-use paper towel.
2. Put the clean and dry parts in a clean storage bag or a clean environment. Do not store parts in an airtight container/bag if moist. It is important that all residual moisture dries.

7. Assembling your Hands-free Collection Cups

Only use clean and dry parts for assembly. Refer to **chapter 6** for cleaning instructions.

7.1 Selecting the correct breast shield size

Make sure to use the correct breast shield size. The breast shield is an essential component for effective pumping, helping to optimize milk flow.

Refer to figure 5 on Chapter 1.

1. 5-01 Select the breast shield size by measuring the diameter of your nipple.
2. 5-02 Select the breast shield size that is slightly larger than your nipple. 21mm (without the O-ring) and 24mm (already applied with the O-ring) size breast shields are included in the pack. 27mm size breast shield is available for purchase. See **chapter 10** for ordering information.

7.2 Assembling the Hands-free Collection Cups

Only use clean and dry parts for assembly. Refer to **chapter 6** for cleaning.

Parts needed:

- Tubing
- Breast shield (with the O-ring)
- Yellow membrane
- Outer shell

CAUTION

To prevent contamination of your milk:

- Wash hands thoroughly with soap and water before touching your breast pump, parts of the Hands-free Collection Cups and breasts.
- Dry your hands with a fresh towel or a single-use paper towel.
- Avoid touching the inside of the Hands-free Collection Cups components.

To ensure proper and safe function of your breast pump:

- Check the Hands-free Collection Cups' components for wear or damage before use. Replace if necessary.
- Always inspect all parts prior to use for cleanliness. If dirty, clean according to **chapter 6**.
- Only use genuine Medela parts. See **chapter 10** for details.

To prevent damage to the breast pump, all components must be completely dry before use.

Refer to figure  6 on Chapter 1.

1.  6-01 Wash hands thoroughly with soap and water before touching your breast pump, parts of the Hands-free Collection Cups and breasts. Dry your hands with a clean towel.
2.  6-02 Confirm that the O-ring is applied to the breast shield.
3.  6-03 Assemble the valve of the yellow membrane on the front of the breast shield.
4.  6-04 Align the disc of the yellow membrane on the designated circular slot using the alignment tab at the top to ensure that it is at the right place. Press around the edges to make sure it is fully assembled.
5.  6-05 Insert the pouring spout into the opening on the shield.
6.  6-06 Push the outer shell and breast shield together tightly, ensuring it is fully sealed all the way around until you hear a snapping sound.
7.  6-07 Connect the flat ends of tubing to your Hands-free Collection Cups.


Information

O-ring fits all breast shield sizes. To disassemble and assemble the O-ring to a different size of breast shield follow the guidance below.

If you would like to use a different size of breast shield, disassemble the already applied O-ring from the 24 mm breast shield starting by pulling out the bottom tab (as shown in the left image below) from the breast shield. To reassemble the O-ring to a different size of breast shield, start from the top (as shown on the right image below) of the breast shield, stretch it over the outer edge of the breast shield. Ensure it goes on with the ridges facing outwards. Press along the groove, checking that it is in place all the way around.

To disassemble:


Pull out the bottom tab from here

To assemble:

Start placing the O-ring from here.

Ensure that the ridges are straight

7.3 Preparing for pumping

Refer to figure 7 on Chapter 1.

1. Assemble the Hands-free Collection Cups as described in **chapter 7.2**.
7-01 Connect the round ends of the tubing to your compatible Medela personal use pump's tubing. For Instructions for Use on how to single or double pump, refer to your Medela compatible personal use pumps instruction for use.

7.4 Placing the Hands-free Collection Cups within your bra

Refer to figure 8 on Chapter 1.

1. 8-01 Ensure that the nipple is centered in the tunnel of the breast shield.
2. 8-02 The Hands-free Collection Cups can be worn under your regular nursing bra. Your bra will hold the Hands-free Collection Cups in place. If you find it uncomfortable or painful during use, try switching to a looser, more stretchy nursing bra.
3. 8-03 The tubing can be connected facing inwards (towards the breasts) or outwards (towards the arms) based on personal preference. Ensure that the tubing is not pinched or bent during use.
4. 8-04 The Hands-free Collection Cups can be used in a vertical or slightly tilted position for better visibility of the nipple and milk ejection. The scale on the Hands-free Collection Cups give you an approximate indication of the volume during pumping.

8. Handling of breast milk

8.1 Removing the Hands-free Collection Cups from your bra

Information

To avoid milk spilling out of the tunnel, when finished, keep the pump running and slightly lean forward to ensure all the milk inside the breast shield tunnel empties into the collection cups.

Refer to figure 9 on Chapter 1.

1. 9-01 While the pump is still running, slightly lean forward and ensure all milk from the nipple tunnel enters the collection container. Next, turn off pump.
2. 9-02 Detach the tubing from the Hands-free Collection Cups.
3. 9-03 Carefully remove the Hands-free Collection Cups from your bra.

8.2 Pouring milk from Hands-free Collection Cups

Refer to figure 9 on Chapter 1.

1. 9-04 Place the Hands-free Collection Cups in an **upright position** on a surface where they will not be easily knocked over.
2. 9-05 **With the rounded outer shell facing down**, pour breast milk through the pouring spout into a bottle or container for storage. **Do not use the Hands-free Collection Cups to store milk. Do not store assembled collection cups in the refrigerator.**
3. Fully disassemble and clean after each use and sanitize once per day as described in **chapter 5** and **chapter 6**.

8.3 General storage guidelines

For information on storing and thawing your breast milk, visit www.breastmilkguidelines.com (U.S.) or <https://www.medela.us/mbus/tips-and-solutions/11/breastmilk-collection-and-storage> (U.S.) <https://www.womenshealth.gov/breastfeeding/pumping-and-storing-breastmilk> (U.S) or <https://www.medela.ca/breastfeeding/moms-journey/storing-and-thawing-breast-milk> (Canada) or <https://www.medela.com/breastfeeding/mums-journey/storing-and-thawing-breast-milk>.

8.4 Long-term storage

Before long-term storage, disassemble and sanitize all parts (except for the tubing). Do not store the Hands-free Collection Cups in assembled state over an extended period of time.

9. Troubleshooting

Problem	Solution
There is low or no suction	<ul style="list-style-type: none">• Make sure that all the Hands-free Collection Cups components are clean and dry and that connections are secure.• Confirm that the O-ring is applied on the breast shield to ensure the suction.• While pumping, make sure the breast shields form a complete seal around the breast.• Refer to the Troubleshooting section in your Medela's personal use pump instruction for use.

10. Ordering information

Accessories

Article number	Product
Depending on your region, check local website/shops for availability.	21mm Breast shields compatible with Medela Hands-free Collection Cups*
	24mm Breast shields compatible with Medela Hands-free Collection Cups*
	27mm Breast shields compatible with Medela Hands-free Collection Cups*
	Membrane compatible with Medela Hands-free Collection Cups
	O-ring compatible with Medela Hands-free Collection Cups
	Tubing compatible with Medela Hands-free Collection Cups and Pump In Style®
	Tubing compatible with Medela Hands-free Collection Cups and Freestyle™ and Swing Maxi™ / Duo

* If the result of expression is unsatisfactory or expression is painful, contact your lactation consultant or breastfeeding specialist. A different size of breast shield may help make expression more comfortable and effective.

i Information

If you have not resolved the problem or in case of difficulties in finding your desired accessory or replacement of lost or defective parts, or any other further questions, please contact Medela customer service.

For contact data, visit www.medela.com, select your country.

For any other Medela products information, visit www.medela.com.

11. Warranty

For US/CA only

This product is warranted by Medela LLC, to the original retail purchaser to be free from defects in material and workmanship for the period of 90 days from the date of purchase. Warranty can only be claimed in the country of purchase. In the event of defect, at Medela's option, Medela will replace this product, without charge for such replacement parts. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

Before you make a claim under this warranty or should you need additional information concerning this warranty, it may save time and expense to call Medela Customer Service (see contacts below). If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, prepaid, together with your dated bill of sale or other proof of purchase and a brief statement of the problem to the following address:

USA
Medela LLC – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS
Medela Customer Service (toll free):
1-800-435-8316

CANADA
Medela Canada
4160 Sladeview Crescent Unit #8
Mississauga, ON L5L 0A1
ATTENTION: RETURNS
info@medela.ca

Contact Customer Service for a return authorization number. Returns are not accepted without an authorization number.

For countries except US/CA:

Find information on the international warranty on www.medela.com/ewarranty.

12. Disposal

Breast shields, outer shells, tubing, o-ring and membranes

The parts are made of plastics that are not harmful to the environment when disposed of as household waste. Recycle or dispose of according to local regulations.

13. Meaning of symbols

The following tables explain the meaning of the symbols found on the product packaging/labeling.

	General safety alert symbol, points to information related to safety. ¹		Identifies the manufacturer. ²
	Indicates the date of manufacturing. ³		Indicates the part number of the product. ⁴
	Read and follow the Instructions for Use. ⁵		Contains fragile goods. Handle with care. ⁶
	Keep away from sunlight. ⁷		Keep away from rain. Keep in dry conditions. ⁸
	Indicates that the package is capable of being recycled. ⁹		Indicates the lot/batch number. ¹⁰

References

- 1 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.2 Symbol 2 General warning sign.
- 2 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.1.1 Manufacturer.
- 3 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.1.3 Manufacturing Date / ISO 7000-2497, Graphical symbols for use on equipment, Date of manufacture.
- 4 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.1.6 Article number / ISO 7000- 2493, Graphical symbols for use on equipment, Catalogue number.
- 5 ISO 15223-1, Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements, Clause 5.4.3, Consult Instructions for Use or consult electronic Instructions for Use / ISO 7000-1641 Operator's manual; operating instructions.
- 6 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.3.1, Fragile, handle with care / ISO 7000-0621, Graphical symbols for use on equipment, Fragile, handle with care.
- 7 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.3.2 Keep away from sunlight / ISO 7000-0624, Graphical symbols for use on equipment, Keep away from sunlight.
- 8 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.3.4, Keep away from rain / ISO 7000-0626, Graphical symbols for use on equipment, Keep away from rain.
- 9 ISO 7000-1135, Graphical symbols for use on equipment, General symbol for recovery/recyclable.
- 10 ISO 15223-1, Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements, Clause 5.1.5 Batch code / ISO 7000-2492, Graphical symbols for use on equipment, Batch code.

14. Technical specifications

Materials touching skin or coming in contact with milk

- Breast shield : Polypropylene
- Outer shell : Polypropylene, Thermoplastic Elastomer
- Membrane : Silicone
- O-ring : Silicone
- Tubing : Silicone, Polycarbonate

Not intentionally made of BPA (Bisphenol A).

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Medela wordmark and logo and Freestyle are registered in the U.S. Patent and Trademark office and the Canadian Intellectual Property Office. Pump in Style is registered in the U.S. Patent and Trademark Office. Science of Care and Swing Maxi are trademarks of Medela.
Medela Canada Inc.