# Instructions for Use

Axiom® Thoracic Catheters and Mediastinal Drains

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## PRECAUTIONS

- The operative site should be dry and free of debris prior to closure.
- Follow universal precautions when inserting and securing the catheters/drains. Special attention is necessary when handling the catheter/drain with instruments. It is not recommended to use wires directly on catheters/drains to hold them in place. When wires are used to realign sternal segments, caution should be taken not to nick or entrap any catheter/drain placed sub-sternally. Proper placement of the catheter/drain in tissue layers and at the exit site should be observed to prevent tube kinking. The catheter/drain may be cut or torn by coming in contact with sharp objects or when subjected to compression or excessive overpressure by a milking roller, etc.
- If catheter/drain is placed on anterior surface of heart after procedure, it should be wiggled or spun before final tightening of sterile wires (or any sutures) to assure that catheter or drain is not caught in wires/sutures and exposed to risk of cutting or compressing when wires/sutures are tightened.
- Fluid retention may result in inefficient evacuation. This could occur as a result of the catheters/drains holes being outside of the tissue layers.
- The eyelets of the catheter/drain must remain inside the thoracic cavity at all times during patient use.
- Surgeon should select the type and size of device that is appropriate for the specific procedure, according to the patient size and effluent volume.
- For devices coated with ClotStop®: For optimum performance, it is recommended that this device should be soaked in sterile water for two (2) minutes prior to use.
- It is common practice to have a catheter/drain exit the wound site via a separate incision – not the main incision – to prevent irritation of incision, possible infection, dehiscence. A surgeon may choose otherwise.
- An adequate number of catheters/drains (at surgeon's discretion) should be used to ensure that all areas will be drained. A separate "Stab Wound" is usually used for each catheter/drain.
- "Stab Wounds" are placed far enough from main incision to allow for repositioning (staged withdrawal) of catheter/drain.
- Catheters/drains usually are combined into one central suction device, but for special reasons, separate suction may be used.

## DIRECTIONS FOR USE

### 1. Handling and Operative Techniques

- Peel open the package and remove the Axiom® catheter/drain using sterile technique. Axiom® catheter/drain coated with ClotStop® include a packaging sheath that must be removed from the catheter/drain prior to patient use.
- Users should be familiar with cardio-thoracic surgical procedures and techniques before using Axiom® thoracic catheters and mediastinal drains. Only qualified healthcare professionals should insert, manipulate and remove thoracic catheters and mediastinal drains.
- The insertion technique, location and catheter/drain size are at the discretion of the physician to optimize drainage efficiency. Care must be taken when inserting the thoracic catheter or mediastinal drain to avoid puncturing the lung and/or other organs.
- The thoracic catheter or mediastinal drain should be secured properly in place to avoid accidental dislodgment or removal.

Made in USA by:  
Axiom® Medical Inc.  
19320 Van Ness Avenue,  
Torrance, CA 90501, USA  
Tel: +1 (800) 221-8569 (USA)  
+1 (310) 533-9020 (Worldwide)  
Fax: +1 (310) 533-8127  
www.AxiomMed.com

SENSS Consulting Ltd  
The Black Church, St Mary's Place  
Dublin 7, D07 P4X, Ireland

Axiom®, Atraum™, ClotStop®, Anchor-Eye®, Stab-Pull®, and Eye-Zone® are registered trademarks of Axiom® Medical, Inc.

© 2020 Axiom® Medical Inc. All rights reserved.

IFU-001-EN Rev A  
Revision date: 2021-04-26

## PRODUCT DESCRIPTION

**Axiom® Thoracic Catheters and Mediastinal Drains** are made from medical grade silicone and are available in a wide variety of sizes (including pediatric) and configurations:

| Round silicone thoracic drains | Flat silicone mediastinal drains |
|--------------------------------|----------------------------------|
| Straight<br>8FR to 40FR        | Right Angle<br>16FR to 40FR      |

- Patient tube connections, insertion site and catheter/drain patency should be checked regularly to confirm proper operation.
- Car should be taken during catheter/drain removal from the patient to avoid damaging the catheter/drain.
- The catheter/drain should be measured after removal and before skin closure to assure that no portion of the device was left in patient due to damage.

## 2. Placement

Once the surgery on the organ and other tissues is complete, the following steps are to be followed:

- An additional "Stab Wound" or "Buttonhole" incision, just slightly larger than the catheter/drain(s), is made a short distance from the main wound incision.
- The catheter(s)/drain(s) are positioned in the wound site where they will work most efficiently with the Stab-Pull® end toward the "Stab Wound".
- The major structures are returned to original positions bringing the catheter(s)/drain(s) around them toward the "Stab Wound".
- A surgical snap or curved clamp is pushed through the "Stab Wound" from exterior to interior of skin thickness.
- The clamp is applied to the Stab-Pull® end and pulls it through the "Stab Wound" from interior to exterior.
- The curved flange on the Stab-Pull® end may be cut off or left in place and the connector tubing is attached to the catheter/drain.
- A suture loops or Axiom® Suture Ring may be used to secure the protruding end of the tube to the patient's skin.
- If applicable, suction is turned on after all junctions are checked for security.

## 3. Removal

- Leaving the soft silicone catheter/drain implanted for any period of time so as to cause tissue around the catheter/drain can interfere with easy removal and affect the performance of the catheter/drain. The surgeon should monitor the patient's rate of wound healing. Gentle flushing of wound site with sterile water should soften encrustation and allow easy, non-irritating, removal of catheter/drain.
- The Axiom® thoracic catheters and mediastinal drains are intended for short term use, i.e. continuous use between 60 minutes and 30 days. Routinely, the devices are used for 24-36 hours, but may be used up to 29 days per physician discretion. The device should be removed per physician's discretion up to 30 days of use.

## DISPOSAL

Dispose of used device per healthcare facility biohazard protocol.

These products are designed for single use only. Discard promptly after single patient use. Do not re-sterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may lead to infection or transmission of blood-borne pathogens to patients and users.

## STORAGE

No special storage conditions required other than to store in a dry, dark, and cool place and to avoid extreme temperatures and humidity. Do not use after expiry date.

## STERILITY

This device has been sterilized, is for single use only, and is not to be reused. As long as the packaging remains sealed and uncompromised the contents within each package are sterile. Do not use after expiration date. Axiom® Medical, Inc. will not be responsible for any products that are resterilized, nor accept for exchange or credit any product that has been opened but not used by the patient or purchaser.

## GLOSSARY OF SYMBOLS

In accordance with:

- ISO 15223-1:2016 (Medical Devices used with – medical device labels, labelling and information to be supplied – Part 1: General requirements)
- ISO 7000:2019 (Graphical symbols for use on equipment)

| Symbol | Symbol Title   | Definition   | Standard Ref. #                                |
|--------|--|--|--|
|        | Manufacturer   | Indicates the medical device manufacturer  | ISO 15223-1 Ref. #5.1.1<br>ISO 7000 Ref. #3082 |
|        | Use-by date  | Indicates the date after which the medical device is not to be used  | ISO 15223-1 Ref. #5.1.4<br>ISO 7000 Ref. #2607 |
|        | Batch code   | Indicates the manufacturer's batch code so that the batch or lot can be identified                               | ISO 15223-1 Ref. #5.1.5<br>ISO 7000 Ref. #2492 |
|        | REF  | Indicates the manufacturer's catalog number that the medical device can be identified.                           | ISO 15223-1 Ref. #5.1.6<br>ISO 7000 Ref. #2493 |
|        | Sterilized using ethylene oxide  | Indicates a medical device that has been sterilized using ethylene oxide   | ISO 15223-1 Ref. #5.2.3<br>ISO 7000 Ref. #2501 |
|        | Sterilized using irradiation   | Indicates a medical device that has been sterilized using irradiation  | ISO 15223-1 Ref. #5.2.4<br>ISO 7000 Ref. #2502 |
|        | Do not re-sterilize  | Indicates a medical device that is not to be re-sterilized   | ISO 15223-1 Ref. #5.2.6<br>ISO 7000 Ref. #2608 |
|        | Do not use if the sterile barrier system or its integrity is compromised | Indicates a medical device that should not be used if the sterile barrier system or its integrity is compromised | ISO 15223-1 Ref. #5.2.8<br>ISO 7000 Ref. #2606 |
|        | Keep away from sunlight  | Indicates a medical device that needs protection from light sources  | ISO 15223-1 Ref. #5.3.2<br>ISO 7000 Ref. #0824 |

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| Symbol | Symbol Title   | Definition   | Standard Ref. #  |
|--------|--|--|--|
|        | Keep dry   | Indicates a medical device that needs to be protected from moisture  | ISO 15223-1 Ref. #5.3.4<br>ISO 7000 Ref. #0626           |
|        | Do not re-use  | Indicates a medical device that is intended for one single use only  | ISO 15223-1 Ref. #5.4.2<br>ISO 7000 Ref. #1051           |
|        | Consult instructions for use   | Indicates the need for the user to consult the instructions for use  | ISO 15223-1 Ref. #5.4.3<br>ISO 7000 Ref. #1641           |
|        | Caution  | To indicate that the current situation needs the user awareness in order to avoid undesirable consequences   | ISO 15223-1 Ref. #5.4.4<br>ISO 7000 Ref. #0434           |
|        | For prescription only  | Caution: Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner.                                       | N/A – Not derived from Standard                          |
|        | Medical device   | Indicates that this item is a medical device   | N/A – From new revision of ISO 15223-1 not yet published |
|        | Unique Device Identification   | Indicates a carrier that contains Unique Device Identifier   | N/A – Not derived from Standard                          |
|        | Quantity   | Indicates the number of units per package  | N/A – Not derived from Standard                          |
|        | European Conformity  | Signifies European CE mark for medical devices. Indicates manufacturer declaration that the product complies with applicable European regulations. | N/A – Not derived from Standard                          |
|        | Authorized representative in the European Community / European Union | Indicates the authorized representative in the European Community / European Union   | ISO 15223-1 Ref. #5.1.2                                  |

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