# Trolley

For Disposable or Reusable Collection System

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HR Upute za upotrebu

# Table of Contents


<b>Overview set-up and accessories</b> .....	<b>3</b>
A Set-up .....	3
B Accessories overview .....	4
C Attaching Disposable Collection System: Serial connection .....	5
D Attaching Disposable Collection System: Change-over valve .....	6
E Attaching Disposable Collection System: Change-over valve and Serial connection .....	7
F Attaching Reusable Collection System: Serial connection .....	8
G Attaching Reusable Collection System: Change-over valve .....	9
<b>1 Warnings and safety instructions</b> .....	<b>10</b>
<b>2 Intended use and Intended user</b> .....	<b>10</b>
<b>3 Description and installation</b> .....	<b>10</b>
<b>4 Before use</b> .....	<b>11</b>
<b>5 Attaching the Fluid Collection Systems</b> .....	<b>12</b>
Disposable Collection System: Serial connection up to 5 jars .....	12
Disposable Collection System: Change-over valve .....	13
Disposable Collection System: Change-over valve and serial connection .....	13
Reusable Collection System: Serial connection up to 5 jars .....	14
Reusable Collection System: Change-over valve .....	14
<b>6 Cleaning and Reprocessing</b> .....	<b>15</b>
<b>7 Warranty</b> .....	<b>16</b>
<b>8 Disposal</b> .....	<b>16</b>
<b>9 Accessories</b> .....	<b>17</b>
<b>10 Technical specifications</b> .....	<b>20</b>
<b>11 Signs and symbols</b> .....	<b>21</b>

## Item overview

077.0188 Trolley A (basic unit 200.8219) with instructions for use 201.0255  
077.0190 Trolley B (basic unit 200.8219) with instructions for use 201.0256

# Overview set-up and accessories

## A Set-up

## B Accessories overview

**Reusable jars**

077.0125: 0.25l  
 077.0155: 0.5l  
 077.0110: 1l  
 077.0120: 2l  
 077.0130: 3l  
 077.0150: 5l

**Reusable lids**

077.0440    077.0420/30

**Reusable sets**

077.0715: 1l  
 077.0716: 2l

**Serial connection**

077.0542

**Change-over valve**

077.0563

**PC jars**

077.0082    077.0085

**Liners**

077.0083/84    077.0086/87  
 101032951

**Serial connection**

077.0160    077.0095

**Vacuum tubing (from suction pump to collection system)**  
**Vacuum tubing, silicone, 7x12 mm**

			077.1019
077.0054: 50cm	077.0921: 30cm		077.1018
077.0901: 100cm	077.0922: 60cm		077.1015/16
077.0902: 200cm	077.0189: 300cm		
077.0905: 500cm			
077.0900: 2500cm	077.0154: 25cm		
	077.0911: 33cm		
077.0931: 25cm	077.0912: 60cm		
077.0185: 60cm	077.0913: 100cm		

**Filters (optional)**

			077.0573	101035263	101035285	101035289
077.0581	077.0092	077.0571	077.0575	101035264	101035286	
			077.0577	101035265	101035287	
			077.0582	101035284	101035288	

## C Attaching Disposable Collection System: Serial connection up to 5 jars

### Necessary parts:

077.0082 (1.5)  
077.0085 (2.5)

077.0083/84 (1.5)  
101032951 (1.5)  
077.0086/87 (2.5)







077.0160/  
077.0158 (100 pcs)

077.0095

077.0189 (3 m)

1. 
2. 
3. 
4. 
5. 
6. 
7. 
8. 
9. 

## D Attaching Disposable Collection System: Change-over valve

### Necessary parts:


077.0082 (1.5)  
077.0085 (2.5)

077.0083/84 (1.5)  
101032951 (1.5)  
077.0086/87 (2.5)

077.0563

077.0189 (3 m)

077.1018

## E Attaching Disposable Collection System: Change-over valve and Serial connection

### Necessary parts:

077.0082 (1.5l)	077.0083/84 (1.5l)	077.0563	077.0160/ 077.0158 (100 pc)	077.0095	077.0189 (3 m)	077.1018
077.0085 (2.5l)	101032951 (1.5l) 077.0086/87 (2.5l)					

- 1.
- 2.
- 3.
- 4.
- 5.
- 6.
- 7.
- 8.
- 9.

# F Attaching Reusable Collection System: Serial connection up to 5 jars

## Necessary parts:

077.0110/20/30/50





077.0420/30/40

077.0542

077.0189 (3 m)

## G Attaching Reusable Collection System: Change-over valve

### Necessary parts:

077.0110/20/30/50

077.0420/30/40


077.0563

077.0189 (3 m)

077.1018

# 1 Warnings and safety instructions

## **WARNINGS**

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

## **WARNINGS**

- For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.
- No modification of this equipment is allowed.

## **CAUTIONS**

Indicates a potentially hazardous situation which, if not avoided, could result in minor injury.

## **CAUTION**

Do not lift the trolley by the handle.

## **Notes**

- Indicating useful information about the safe use of the device.
- The trolley is approved exclusively for use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when the trolley is used in combination with the original Medela accessories.
- These instructions for use contain essential information on using this equipment safely and effectively. Before use, thoroughly review these instructions for use and the manuals of all equipment which will be used during the procedure and use the equipment as instructed. Keep these instructions for use with the trolley for later reference.
- If you need further assistance with device operation, please contact your local Medela customer service representative.

# 2 Intended use and Intended user

## **Intended use/indications**

The Medela Trolley is intended to be used throughout the hospital to hold up to 5 suction jars (PSU or PC with disposable liners). It is designed for indoor use and equipped with 5 wheels. Do not use the equipment for any purpose other than its intended application.

## **Intended user**

The trolley should only be operated by properly trained staff.

**CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.**

**Any serious incident that has occurred in relation to the device should be reported to Medela AG and the relevant Competent Authority.**

## 3 Description and installation

See pictures chapter A Set-up, page 3.

**On initial delivery check the delivery package of the trolley for completeness and that the parts are intact**

1. 5 castors
2. star base
3. trolley pole
4. jar holder
5. handle

**Assembling steps**

1. Attach castors to star base.
2. Attach pole. Screw thread is facing up.
3. Attach jar holder.
4. Secure jar holder with handle.

**Activating/deactivating the brakes**

To activate, press lever down.

To deactivate, press lever up.

## 4 Before use

**Perform the following safety checks**

- Examine the trolley for obvious damage or safety defects before use.
- Prior to any use examine all accessories and replace if necessary:
  - suction jars, lids change over valves and liners for cracks, brittle and flawed spots
  - tubing for cracks and brittle areas
- Make sure that connectors are firmly attached.

As a system check, evacuate the complete system (including jars) at maximum vacuum before actual use. This enables you to confirm the set-up is tight.

**CAUTION**

The trolley must remain in an upright position during use.

# 5 Attaching the Fluid Collection Systems

Various options are available to set up the Disposable and Reusable Collection System.

**WARNING:** For the disposable liner system it is important to close all ports with the attached caps before removing vacuum.

## 5.1 Disposable Collection System: Serial connection up to 5 jars (chapter C)

**See pictures chapter C, page 5.**

1. Insert the number of desired jars including liners into the jar holder.
2. Cut off 10 cm of the tubing of the serial vacuum connection.
3. Insert serial vacuum connections into the suction port of all jars except for the right lower jar.
4. Insert coupling piece into the right lower (= last) jar.
5. Connect all serial vacuum connections to the next suction port (connect them to the vacant Y-piece of the next serial connection or the coupling piece of the last jar).
6. a. Pull the suction tubing from the vacuum source from below the jar holder through the hole closest to jar and b. connect to the suction port (Y- piece) of the lower jar. Connect this suction tubing to suction source. We recommend using the Medela Dominant Flex suction pump.
7. Insert serial patient connections into all serial ports except for the right lower (=last) jar.
8. Insert coloured angle piece to the patient port of all liners. Connect open end of the serial patient connections with the coloured angle piece of the next jar's patient port.  
(Serial Port Jar 1 to patient port Jar 2 ...)
9. Attach desired patient tubing to angle piece of lower jar.

## 5.2 Disposable Collection System: Change-over valve (chapter D)

**See pictures chapter D, page 6.**

### **Note**

When using change-over valve instead of serial connection, suction is switched from one jar to the other manually. No reduction in flow rate is to be expected.

1. Insert a jar including a liner to each side of the change-over valve.
2. Connect the left suction tubing from the change-over valve to the suction port of the left jar.  
Do the same for the right side.

3. a. Take suction tubing from the vacuum source and pull from below the jar holder through the hole closest to the change-over valve and connect it to the middle port of the change-over valve.
- b. Connect the opposite end of this suction tubing to vacuum source. We recommend using the Medela Dominant Flex suction pump.
4. If necessary, insert coloured angle piece to the patient port(s). Connect desired patient tubing to the left and the right jars.

## 5.3 Disposable Collection System: Change-over valve and serial connection (E)

**See pictures chapter E, page 7.**

### **Note**

When using serial connection, the flow rate of the suction source is reduced. Make sure the flow rate is still adequate for your desired application.

1. Insert 4 jars including liners to each side of the change-over valve.
2. Cut off 10cm of the tubing of the serial vacuum connection.
3. Insert two serial suction connections to the jars next to the change-over-valve with its open tubing facing the upper jars.
4. Remove the two coupling pieces of both the change-over-valve tubing and plug these coupling pieces into the suction port of the two upper jars.
5. a. Connect the suction tubing from the change-over-valve with the Y piece in the vacuum port of the neighboring jars.
- b. Also connect the serial vacuum connections from the middle jars to the upper jars.
6. a. Take the suction tubing from the vacuum source and pull from below the jar holder through the hole closest to the change-over valve.
- b. Connect the suction tubing to the middle port of the change-over-valve.  
    Connect the suction tubing to suction source. As a suction source we recommend using the Medela Dominant Flex suction pump.
7. Insert two serial patient connections into serial port of left lower jar and right upper jar.
8. Insert coloured angle piece into the patient ports of the liners. Connect the serial patient connections with the neighboring jars.
9. Attach desired patient tubing to coloured angle piece of first and last jars.

## 5.4 Reusable Collection System: Serial connection up to 5 jars (chapter F)

**See pictures chapter F, page 8.**

### **Note**

When using serial connection, the flow rate of the suction source is reduced. Make sure the flow rate is still adequate for your desired application.

1. Insert required number of jars into jar holder. Remove mechanical overflow protection from all jars except the front jar.
2. Attach lids and put a reusable serial connector in the vacuum ports of all jars except the front jar. The tubing shall point towards the previous jar.
3. Connect the tubing of the reusable serial connectors with the patient port of the previous jar.
4. a. Take suction tubing from the vacuum source and pull from below the jar holder through the hole closest to the front jar.  
b. Plug in suction tubing into vacuum port on lid of front jar.  
c. Connect desired patient tubing to the patient port on the right lower (= last) jar.

## 5.5 Reusable Collection System: Change-over valve (chapter G)

**See pictures chapter G, page 9.**

### **Note**

When using change-over valve, suction is switched from one jar to the other manually. No reduction in flow rate is to be expected.

1. Insert a jar, having the overflow protection attached, to each side of the change-over valve.
2. a. Connect the left tubing from the change-over valve to the left jar. Do the same for the right side.  
b. Connect suction tubing from the vacuum source to the change-over valve, leading the suction tubing from below through a hole in the jar holder to the jar.
3. Attach desired patient tubing to the patient ports of the left and right jars.

## 6 Cleaning and Reprocessing

### Cleaning and Reprocessing Instructions

Reusable devices and components from Medela are delivered non-sterile and are intended for reuse. Prior to use clean and reprocess the product following the procedure described below.

The following cleaning and reprocessing processes have been validated according to international standards and guidelines:

- ISO 17664:2004 (EN) Sterilization of medical devices — Information to be provided by the manufacturer for the processing of resterilizable medical devices
- Guidance for Industry and FDA Staff: March 17, 2015 Processing/Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling

It is the responsibility of the user/processor to ensure that the cleaning/reprocessing procedure is performed using equipment, materials and personnel which are suitable to ensure the effectiveness of the process. Any deviation from the following instructions should be validated by the user/processor to ensure effectiveness of the process.

Materials and equipment:

- Personal Protective Equipment (PPE: disposable gloves, proper protective gear)
- Lint free nonwoven wipes

*Initial cleaning:* Perform initial cleaning directly after the use of the device (within a maximum of 2 hours after procedure). Wipe the surface of the device to remove all gross debris by using a lint free nonwoven wipe moistened with water (< 40°C, < 104°F).

*Thorough manual cleaning:* Use a lint free nonwoven wipe moistened with tap flowing water (< 40°C, < 104°F).

### Note

Never use steel brushes or steel wool for cleaning.

Water temperatures higher than 45°C (113°F) can cause protein coagulation which may negatively affect the efficacy of the cleaning process.

*Manual disinfection:* Use only disinfectant wipes with active substances such as alcohols, quaternary ammonium compounds (Quats), hypochlorites, hydrogen peroxides, glutaraldehydes. For effective and safe disinfection please follow the instructions of the disinfectant manufacturer.

*Drying:* After manual cleaning and disinfection, store the product on a clean and dry cloth and let it dry or follow instructions for use of the disinfectant manufacturer.

*Inspection:* After reprocessing according to the steps described above, the device and its components must be inspected for signs of degradation that may limit the useful life and / or performance of the device, such as the following: visible corrosion, mechanical wear, abrasion, damage, or deformation. Discard the device if any signs of degradation are evident.

*Storage:* Store the product dry and dust-free.

## 7 Warranty

Medela AG guarantees the trolley to be free of defects in materials and workmanship for a period of 1 year from the date of delivery. Should any material defects appear, the appropriate parts will be replaced free of charge during this period, provided it does not result from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with the Medela trolley.

In no event shall Medela AG be liable for claims which exceed the scope of the warranty described, including liability for consequential damage, etc. The right to the replacement of faulty parts will not be recognised by Medela if any work has been carried out on the trolley by unauthorised persons. This warranty is subject to the appliance being returned to a Medela service centre.

## 8 Disposal

The trolley is made of metals and plastics and should be disposed of in accordance with the European directives 2002/95/EC and 2002/96/EC. In addition, follow your local disposal guidelines when disposing of the trolley and accessories.

## 9 Accessories

<b>Filters with overflow protection</b>	Bacteria filter with overflow protection	
	077.0571	Medela connections
	077.0572	Conical connections
	077.0577	Medela and conical connection
	Virus filter with overflow protection	
	101035263	Medela connections
	101035264	Conical connections
	Bacteria/odour filter with overflow protection	
	077.0092	Medela connections
	101035285	Conical connections
	Virus/odour filter with overflow protection	
	101035286	Medela connections
101035287	Conical connections	
<b>Filters without overflow protection</b>	Bacteria filter	
	077.0573	Medela connections
	077.0574	Conical connections
	Virus filter	
	101035265	Medela connections
	101035284	Conical connections
	Bacteria/odour filter	
	077.0575	Medela connections
	077.0576	Conical connections
	Virus/odour filter	
	101035288	Medela connections
	101035289	Conical connections
	Particle filter	
	077.0582	Medela connections
	<b>Filters</b>	077.0581
077.0582		Disposable particle filter with Medela connections
<b>Reusable sets with PSU suction jar and lid</b>	077.0715	1 l Set
	077.0716	2 l Set
<b>Disposable suction liners</b>	077.0083	1.5 l
	077.0084	1.5 l with solidifier
	101032951	1.5 l with solidifier
	077.0086	2.5 l
	077.0087	2.5 l with solidifier
<b>Suction jars PC (for disposable suction liners)</b>	077.0082	1.5 l
	077.0085	2.5 l

<b>Silicone suction tubing Ø 7 x 12 mm (leading to suction source)</b>	077.0902	w/o coupling pieces, 2 m
	077.0912	with 1 coupling piece, 60 cm
	077.0900	w/o coupling pieces, 25 m
	077.0901	w/o coupling pieces 1 m
	077.0905	w/o coupling pieces, 5 m
	077.0913	with 1 coupling piece, 100 cm
	077.0921	with 2 coupling pieces, 30 cm
	077.0189	with 2 coupling pieces, 3 m
	077.0922	with 2 coupling pieces, 60 cm
077.0931	with 1 coupling piece and 1 coupling piece 90°	
<b>Angle pieces: attach to lid of Reusable Collection System (patient port)</b>	077.1015	Ø 6–10 mm
	077.1016	Ø 10–14 mm
<b>Coupling piece: Used to connect tubing to vacuum port of pump</b>	077.1018	Coupling piece
	077.1019	Coupling piece 90°
<b>Patient tubing</b>	077.0951	Disposable (PVC) tubing, 180 cm, with fingertip, sterile
	077.0963	PVC tubing Ø 12.7 x 18 mm, 200 cm for patient connection Ø 10–14 mm
	077.0965	Disposable PVC tubing, 150 cm CH 25, sterile
	077.0970	Silicone tubing Ø 5 x 10 mm, 25 m
<b>Reusable suction jars, polysulfone</b>	077.0110	1 l
	077.0120	2 l
	077.0130	3 l
	077.0150	5 l
	077.0125	0.25 l
	077.0155	0.5 l
<b>Reusable lids for suction jars, polysulfone, with conical patient connection and overflow protection</b>	077.0440	Small lid, Ø 6–10 mm for 0.25 l and 0.5 l suction jars
	077.0420	Large lid Ø 6–10 mm
	077.0430	Large lid Ø 10–14 mm
<b>Two or more jars (Disposable Collection System)</b>	077.0095	Serial connection, vacuum
	077.0160	Serial connection, patient, single
	077.0158	Serial connection, patient, (100 pcs)
	077.0563	Change over valve with clampholder and tubing set
<b>Two or more jars (Reusable Collection System)</b>	077.0563	Change over valve with clampholder and tubing set
	077.0542	Reusable serial connector to connect max. five 5 l jars

**WARNING**

If the trolley is used with third party patient interfacing devices (e.g. cannulas, catheters), they must:

- have a CE mark and (if necessary) local registration
- be compatible with Medela accessories without negatively impacting safety and the performance of the trolley

Tubing connection on reusable lids of jars: Ø 6–10 mm, Ø 10–14 mm

Tubing connection on disposable liners: Ø 7–11 mm

**WARNING**

When combining Medela parts and a third party interfacing device, the user takes responsibility for the entire system and should test the combination to ensure that the performance and safety are properly maintained.


# 10 Technical specifications


Weight: 4 kg/8.8 lbs  
(without packaging)

HxWxD  
850x600x600 mm

Storage

Operation

# 11 Signs and symbols

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates to follow instructions for use.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates the weight of the device.

This symbol indicates the temperature limitation for operation and storage.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device dry.

This symbol indicates the humidity limitation for operation and storage.

This symbol indicates the atmospheric pressure limitation for operation and storage.

This symbol indicates do not use the device if package is damaged.

This symbol indicates the dimensions of the system.

This symbol indicates to not dispose of the device together with unsorted municipal waste (for EU only).

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates a carton package.

This symbol indicates to keep the device away from sunlight.

This symbol indicates the manufacturer.

This symbol indicates manufacturer's catalogue number.

This symbol indicates manufacturer's batch number.

This symbol indicates the item is a medical device.

This symbol indicates a prescription device. U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).

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# Sadržaj

<b>Overview set-up and accessories</b> .....	<b>3</b>
A Set-up .....	3
B Accessories overview .....	4
C Attaching Disposable Collection System: Serial connection .....	5
D Attaching Disposable Collection System: Change-over valve .....	6
E Attaching Disposable Collection System: Change-over valve and Serial connection .....	7
F Attaching Reusable Collection System: Serial connection .....	8
G Attaching Reusable Collection System: Change-over valve .....	9
<b>1 Upozorenja i sigurnosne upute</b> .....	<b>143</b>
<b>2 Namjena i predviđeni korisnik</b> .....	<b>143</b>
<b>3 Opis i instalacija</b> .....	<b>144</b>
<b>4 Prije upotrebe</b> .....	<b>144</b>
<b>5 Pričvršćivanje sustava za prikupljanje tekućine</b> .....	<b>145</b>
Jednokratni sustav za prikupljanje: Serijsko priključivanje do 5 posuda (poglavlje C) .....	145
Jednokratni sustav za prikupljanje: Preklopni ventil (poglavlje D) .....	146
Jednokratni sustav za prikupljanje: Preklopni ventil i serijsko priključivanje (poglavlje E) .....	146
Višekratni sustav za prikupljanje: Serijsko priključivanje do 5 posuda (poglavlje F) .....	147
Višekratni sustav za prikupljanje: Preklopni ventil (poglavlje G) .....	147
<b>6 Čišćenje i ponovna obrada</b> .....	<b>148</b>
<b>7 Jamstvo</b> .....	<b>149</b>
<b>8 Odlaganje</b> .....	<b>149</b>
<b>9 Pribor</b> .....	<b>150</b>
<b>10 Tehničke specifikacije</b> .....	<b>152</b>
<b>11 Znakovi i simboli</b> .....	<b>153</b>

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Medela UK Ltd.
Medela LLC