



# Invia<sup>®</sup> Abdominal Dressing Kit

## NEGATIVE PRESSURE WOUND THERAPY

- EN Instructions for use
- ES Instrucciones de uso
- FR Mode d'emploi
- DE Gebrauchsanweisung
- IT Istruzioni per l'uso

# Table of contents

- Device description ..... 4
- Intended use and therapy safety information ..... 4
- Indications for use ..... 5
- Contraindications ..... 5
- Warnings and safety instructions ..... 5
- Dressing application ..... 10
- Dressing changes..... 13
- Dressing removal..... 13
- Signs and symbols ..... 14

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician.

**Only for use with Invia Negative Pressure Wound Therapy Systems**

## Device description

The Invia Abdominal Dressing Kit consists of an Organ Contact Layer (OCL), two foam pads, four transparent films and a suction interface (Invia FitPad).

- The Organ Contact Layer (OCL) is an oval polyurethane film with fenestrations intended to protect the abdominal content and enable fluid drainage from the abdomen.
- The foam is an oval hydrophobic polyurethane foam with open cell structure intended to be placed over the OCL in order to distribute the negative pressure across the wound surface and allow passage of fluids and exudates through to the negative pressure system.
- The transparent film is a thin transparent polyurethane film with acrylic adhesive intended to fixate the wound filler and seal tight to the skin.
- The Invia FitPad is intended to distribute negative pressure to the wound and transport exudates from the abdominal cavity to the canister.

(201.0697) Kit contents		
Description	QTY in kit	Dimensions
Organ Contact Layer (OCL)	1	oval (60x80 cm)
Foam dimensions (Length x Width x Thickness)	2	38x25x2 cm
Transparent film	4	32x26 cm
Suction pad (FitPad)	1	-

### Ordering information:

Invia Abdominal Dressing Kit

**REF** 087.6250 (5 pcs of 201.0697)

## Intended use and therapy safety information

The Invia Abdominal Dressing Kit is intended for use in conjunction with the Invia Liberty Negative Pressure Wound Therapy (NPWT) system. Users are directed to the Invia Liberty NPWT system labeling for additional safety information and instruction for use. To help ensure safe and effective use, the Invia Abdominal Dressing Kit is to be used only with the approved therapy units.

The components of the Invia Abdominal Dressing Kit are packaged sterile and are for single use only. Do not use if sterile package is damaged or opened prior to use. To apply Invia Abdominal Dressing Kit, use clean/aseptic or sterile techniques in accordance with local protocol.

The intended use of the Invia Abdominal Dressing Kit together with Invia Liberty Negative Pressure Wound Therapy system is to manage adult patients with open abdominal wounds with exposed viscera and organs, and including but not limited to adult patients with abdominal compartment syndrome. Do not use the Invia Abdominal Dressing Kit on neonates, children, nursing and or pregnant women since they are excluded from the intended patient population.

It is intended for use in acute hospital settings (trauma, general and plastic surgery wards) and should ideally be applied in the operating theatre by a trained physician.

## Indications for use

The Invia Abdominal Dressing Kit is indicated for temporary bridging of abdominal wall openings where primary closure is not possible and/or repeat abdominal entries may be required. Its intended use is with patients who have open abdominal wounds with exposed viscera and organs, and including but not limited to patients with abdominal compartment syndrome. It is intended for use in acute hospital settings (trauma, general and plastic surgery wards) and should ideally be applied in the operating theatre. The dressing kit is intended for use together with the Invia Liberty NPWT System.

## Contraindications

- Direct positioning of NPWT foam over exposed organs, large veins and arteries, anastomotic sites, tendons or nerves (unless suitably covered with the Organ Contact Layer)
- Necrotic tissue with eschar present
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Malignancy in wound

## Warnings and safety instructions

### WARNING

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

### CAUTION

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

### Safety related tip

Indicating useful information about the safe use of the device.

### WARNINGS

**Invia Abdominal Dressing Kit is intended to be used in conjunction with the Invia Liberty Negative Pressure Wound Therapy system.**

- Please refer to the Invia Liberty Negative Pressure Wound Therapy system product instructions for use prior to initiating therapy.

### Use of Organ Contact Layer (OCL):

- Use of the OCL is necessary to protect exposed bowels from the abdominal wall and from contact to the foam. Preventing adhesions or obstructions that may otherwise form during open abdominal wound management is a critical parameter in achieving timely primary facial closure and reducing the chance of fistula development.

- Make sure the OCL is not cut too small.
- The OCL must completely cover all exposed viscera, prior to application of wound filler and subsequent activation of the Invia Liberty NPWT system.

**Foam placement:**

- Do not tear or cut the foam above the abdominal cavity. Remove any fragments from the foam edges to protect the abdominal cavity from loose foam particles.
- Insufficient filling may cause leakage of fluids onto surrounding skin. Use additional foam if necessary.
- Make sure that the foam does not come into contact with viscera or surrounding skin
- Do not place any foam dressing pieces into blind/unexplored tunnels.
- Always count the total number of pieces of foam used in the wound and document that number on the transparent film and in the patient's chart.
- Fixation of the foam (sutures, staples etc) may cause foam particles to enter the abdominal cavity and/or the margins of the abdominal cavity.

**Film application:**

- Stretching the film during application may cause the film to loosen when negative pressure is applied.

**Intra-abdominal pressure monitoring:**

- Frequent, standardized measurements of intra-abdominal pressure (IAP) and/or abdominal perfusion pressure (APP) are recommended before, during and after treatment with Invia abdominal NPWT. If intra-abdominal hypertension (IAH) is observed or suspected, turn off power to the NPWT system. After full expansion of the foam, obtain a new intra-abdominal pressure measurement.
- If IAH persists, discontinue the use of NPWT and address the underlying condition.

**Bleeding:**

- With or without using therapy, certain patients are at high risk of bleeding complications. **Patients with abdominal wounds must be closely monitored for bleeding as these wounds may contain hidden blood vessels which may not be readily apparent.**
- **Hemostasis must be achieved prior to dressing placement.**

The following types of patients are at an increased risk of bleeding, which, if uncontrolled, could be potentially fatal:

- Patients who have weakened or friable blood vessels or visceral organs in the abdominal wound as a result of, but not limited to:
  - Suturing of the blood vessel (native anastomosis or grafts)/organ
  - Infection
  - Trauma
  - Radiation
- Patients without adequate wound hemostasis
- Patients who have been administered anticoagulants or platelet aggregation inhibitors
- Patients who do not have adequate tissue coverage over vascular structures

- If active bleeding develops suddenly or in large amounts during therapy, or if frank (bright red) blood is seen in the tubing or in the canister, immediately stop therapy, leave dressing in place, take measures to stop the bleeding and seek immediate medical assistance.
- The therapy units and dressings should not be used to prevent, minimize or stop vascular bleeding.
- Do not resume the use of the therapy system until adequate hemostasis has been achieved and the patient is not at risk of continued bleeding.

#### **Protect vessels and organs:**

- All exposed or superficial vessels and organs in the abdominal wound must be completely covered and protected prior to the administration of therapy. **Always ensure that foam dressings do not come in direct contact with vessels or organs.**

#### **Infected blood vessels:**

- Infection may erode blood vessels and weaken the vascular wall which may increase susceptibility to vessel damage through abrasion or manipulation. **Infected blood vessels are at risk of complications, including bleeding, which, if uncontrolled, could be potentially fatal. Extreme caution should be used when therapy is applied in close proximity to infected or potentially infected blood vessels.**

#### **Hemostasis, anticoagulants and platelet aggregation inhibitors:**

- Patients without adequate wound hemostasis have an increased risk of bleeding, which, if uncontrolled, could be potentially fatal. These patients should be treated and monitored in a care setting deemed appropriate by the treating physician.

#### **Hemostatic agents applied at the wound site:**

- Non-sutured hemostatic agents (for example, bone wax, absorbable gelatin sponge, or spray wound sealant) may, if disrupted, increase the risk of bleeding, which, if uncontrolled, could be potentially fatal. Protect against dislodging such agents. Consideration should be given to the negative pressure setting and therapy mode used when initiating therapy.

#### **Sharp edges:**

- Bone fragments or sharp edges could puncture protective barriers, vessels, or organs causing injury. Any injury could cause bleeding, which, if uncontrolled, could be potentially fatal.
- Sharp edges or bone fragments must be eliminated from the wound area or covered to prevent them from puncturing blood vessels or organs before the application of therapy.

#### **Infected wounds:**

- Infected abdominal wounds should be monitored closely and may require more frequent dressing changes than non-infected wounds, dependent upon factors such as wound conditions and treatment goals.
- Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge).

Non-observance can lead to considerable danger to the patient.

**Osteomyelitis:**

- The therapy system should NOT be initiated on a wound with untreated osteomyelitis. Consideration should be given to thorough debridement of all necrotic, nonviable tissue, including infected bone (if necessary), and appropriate antibiotic therapy.

**Dressings removal:**

- Abdominal dressing kit components are not bioabsorbable.
- **Always remove all dressing components from the abdomen at every dressing change.**
- **Count the total number of foam pieces removed from the abdominal wound and ensure the same number of foam pieces was removed as placed.**

**Acrylic adhesive:**

- The Invia Transparent Film has an acrylic adhesive coating, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives.
- If a patient has a known allergy or hypersensitivity to such adhesives, do not use the therapy system.
- If any signs of allergic reaction or hypersensitivity develop, such as redness, swelling, rash, urticaria, or significant pruritus, discontinue use and consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, seek immediate medical assistance.

**Defibrillation:**

- Remove the dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and/or patient resuscitation.

**Magnetic Resonance Imaging (MRI):**

- Therapy Unit: Invia NPWT systems are not for use in the Magnetic Resonance (MR) Environment, so do not take an Invia NPWT device into this environment.
- Dressings: After disconnecting from the Invia NPWT system, dressings can typically remain on the patient with no risk in a MR environment (all components are electrically nonconductive and nonmagnetic items).

**Patients with Spinal Cord Injury** (sympathetic nervous system stimulation):

- In the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation in response to stimulation of the sympathetic nervous system), discontinue Invia NPWT therapy to help minimize sensory stimulation and seek immediate medical assistance.

**Course of therapy:**

- Invia Negative Pressure Wound Therapy Instructions advise 24 hours therapy without interruption.
- If therapy is discontinued for more than 2 hours, the dressing should be replaced and therapy restarted by a physician. If the therapy needs to be interrupted, the tubing should be clamped and the ends of the tubing protected.

 **CAUTIONS**

The following statements describe conditions that may require special care for the safe and effective use of the Invia Abdominal Dressing Kit.

- Patients untreated for malnutrition.
- Noncompliant or combative patients.
- Wounds that involve an enteric fistula.
- To minimize the risk of bradycardia, NPWT must not be placed in the proximity to the vagus nerve.
- Avoid circumferential dressing application.
- The device has not been studied on pediatric patients.
- Patient size and weight should be considered when prescribing this device. Initial lower negative pressure should be considered for certain small or elderly patients who are at risk of fluid depletion or dehydration.
- As a condition of use, this device should only be used by qualified and authorized personnel. The user must have the necessary knowledge of the specific medical application for which Invia Abdominal Dressing Kit is being used.
- When treating abdominal wounds, the therapy should remain on in continuous mode for the duration of the treatment.
- Ensure that tubing is installed completely and without any kinks to avoid leaks or blockages in the suction circuit.
- Maintain regular monitoring of the wound site during therapy to ensure therapeutic treatment and patient comfort.
- As with all adhesive products apply and remove the dressing carefully from sensitive or fragile skin to avoid skin tripping, especially after frequent dressing changes.
- Do not use if packaging is breached or damaged.


# Dressing application

## WARNING


Review all Invia Abdominal Dressing Kit safety information before beginning with the preparation of open abdominal wound. Ensure sufficient hemostasis has been achieved prior to dressing application (refer to Bleeding section under WARNINGS).

## Step 1: Prepare the open abdominal wound

- Eliminate any sharp edges or bone fragments from wound area.
- Ensure any areas of necrosis are appropriately debrided.
- Irrigate abdominal wound as needed.
- Clean and thoroughly dry the periwound area.

## Step 2: Apply the Organ Contact Layer (OCL)


- Remove kit contents from pouch and prepare the OCL in a sterile field.
- If necessary, cut the OCL to a size that fits the abdominal cavity.
- If cutting the OCL to a different size, ensure that each piece removed has been disposed of properly, away from the open wound.
- Center the OCL over the abdominal cavity. Gently work it down to cover all exposed organs, including intestines and liver, and into the paracolic gutters. If there is excess material of the OCL this can be folded back onto itself.
- Make sure that all viscera are protected by the OCL.



 Either side of the OCL may be applied to exposed organs.

 Moistened gloves (saline solution or sterile water) may facilitate application of the OCL.

### Step 3: Apply the foam


- Tear or cut the foam into an appropriate size corresponding to the dimensions of the abdominal cavity.


- Carefully place the foam in the abdominal cavity.
- Sufficiently fill the cavity without over-packing.

### Step 4: Apply the transparent film


Clean and thoroughly dry the periwound area. Consider use of skin preparation product to protect the periwound skin from exposure to moisture and adhesive.


a) Peel the central layer "1"

b) Place the adhesive face down over the foam and remove side layers "1"

c) Peel top layer "2"


d) Detach handling bars along perforated lines


More than one piece of film can be used to cover the wound.

- ☞ When using multiple pieces of film ensure the edges overlap by a minimum of 5–7 cm.
- ☞ Avoid stretching or pulling the film to minimise tension or trauma to the peri-wound skin.
- ☞ Film should extend 5–7 cm beyond wound margin to facilitate adequate seal.
- ☞ Film should be securely anchored to periwound area to maintain an air tight seal.
- ☞ Film can be cut into multiple pieces for easier handling. Make sure to retain a portion of the handling bars on each piece.

### Step 5: Attach the suction pad (FitPad)

- Select an appropriate location over the dressed wound bed for the suction pad (FitPad) to be applied.
- Pinch the transparent film over the area, lifting slightly.
- Cut a small hole in the pinched area (approx. 1 cm diameter).



- Peel off the backing of the suction pad.
- Center the suction pad over the previously cut hole in the dressed wound and apply by pressing firmly for adherence.
- Pay attention to tubing positioning to allow for optimal flow: avoid placement over bony prominences.


### Step 6: Connect the dressing tubing to the pump tubing

 Make sure both parts of the Quick-connector are aligned correctly as shown in the picture below.

 Push the Quick-connector together until you hear a click.

### Step 7: Turn on the pump

- Turn on the pump and select the prescribed therapy setting.
- Finished dressing should be firm to the touch.

 In order to avoid an early trigger of the “system leak” alarm, the high leakage mode can be selected on the menu of the Invia Liberty NPWT pump.

 The pressure level shall always be set according to clinician’s instruction. For guidance the recommended pressure level is  $-125\text{mmHg}$  and constant mode.

## Dressing changes

- Dressings Organ Contact Layer (OCL), foam, transparent film and FitPad) should be changed no later than 48–72 hours.
- The frequency of dressing changes should be based on an evaluation of patient condition rather than standard recommendations.
- In the event of heavy drainage, drainage with sediment or infected wounds, more frequent dressing changes may be needed.
- Check dressings regularly and monitor the wound to check for signs of infection. If there are any signs of systemic infection or advancing infection at the wound site, contact the treating physician immediately.

 The total therapy time should not exceed 30 days.

## Dressing removal

- a. Turn the pump off.
- b. Close the clamp on the dressing tubing.
- c. Close the clamp on the pump tubing.
- d. Disconnect pump tubing from dressing tubing by pressing on the sides of the Quick-connector.
- e. Remove the transparent film from the skin by gently pulling parallel with the skin.
- f. Remove the foam from the wound.

**Ensure the same number of foam pieces that were placed in the wound has been removed.**

 If foam adheres to wound margins, apply normal saline into the dressing and wait for 15–30 minutes before gently removing the foam.

g. Carefully lift the Organ Contact Layer (OCL) out of the abdominal cavity.

**Check that the OCL is intact and that no pieces are left in the abdominal cavity**

- h. Dispose of the dressing components in accordance with local guidelines.
- i. Prepare the wound for the next dressing as described under “Prepare the open abdominal wound”.

# Signs and symbols

This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).

This symbol indicates the temperature limitation for operation and storage.

This symbol indicates the manufacturer.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates manufacturer's batch code.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the date of manufacture.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates to follow instructions for use.

This symbol indicates the humidity limitation for operation and storage.

This symbol indicates do not use the device if package is damaged.

This symbol indicates manufacturer's catalogue number.

This symbol indicates safety related tip.

This symbol indicates the atmospheric pressure limitation for operation and storage.

This symbol indicates number of items.

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Medela AG
Medela UK Ltd.
Medela LLC