



# Medela NICU Solutions


TwistLok® Connectors  
Product Ordering Guide


neonatal solutions




# Ordering Information

Enteral Syringes		Description	Item#	Quantity
	TwistLok 5 mL Syringe	NTL05	100 per case	
	TwistLok 10 mL Syringe	NTL10	100 per case	
	TwistLok 30 mL Syringe	NTL30	50 per case	
	TwistLok 60 mL Syringe	NTL60	50 per case	

Enteral Syringes (bonded TwistLok® connector)		Description	Item#	Quantity
	Intopo (Slip Fit) 1 mL Syringe	INT01	100 per case	
	TwistLok 3 mL Syringe	TLINT03	100 per case	
	TwistLok 5 mL Syringe	TLINT05	100 per case	
	TwistLok 10 mL Syringe	TLINT10	100 per case	
	TwistLok 20 mL Syringe	TLINT20	100 per case	
	TwistLok 30 mL Syringe	TLINT30	50 per case	
	TwistLok 60 mL Syringe	TLINT60	50 per case	


Enteral Extension Sets TwistLok Female Connector, Orange-Striped Minibore Tubing, Non-Removable Slide Clamp (White), Enteral Feeding Direction Tag, TwistLok Male Connector		Length/Priming Vol.	Item#	Quantity
	60", 1.3 mL	TL06018	100 per case	
	60", 2.0 mL	TL060182	100 per case	
	70", 1.50 mL	TL07018	100 per case	

Enteral Extension Sets with Medication Port Two TwistLok Female Connectors, Orange-Striped Low-Volume Minibore Tubing, Two Non-Removable Slide Clamps (White), Medication Port, Non-Vented Plug, TwistLok Male Connector		Length/Priming Vol.	Item#	Quantity
	60", 1.7 mL	TL06019	100 per case	
	60", 2.2 mL	TL060192	100 per case	
	70", 1.90 mL	TL07019	100 per case	

# Ordering Information

## PVC Enteral Feeding Tubes

(For use up to 3 days)



Length/Priming Vol.	Size	Item#	Quantity
16", 0.30 mL	5 French	ENPV1650	50 per case
26", 0.45 mL	5 French	ENPV2650	50 per case
36", 0.90 mL	5 French	ENPV3650	50 per case
16", 0.60 mL	6.5 French	ENPV1665	50 per case
26", 1.00 mL	6.5 French	ENPV2665	50 per case
36", 1.20 mL	6.5 French	ENPV3665	50 per case
16", 1.00 mL	8 French	ENPV1680	50 per case
26", 1.40 mL	8 French	ENPV2680	50 per case
42", 2.40 mL	8 French	ENPV4280	50 per case

## Silicone Enteral Feeding Tubes

(For use up to 30 days)



Length/Priming Vol.	Size	Item#	Quantity
40cm (15.7"), 0.45 mL	5 French	ENSI4050	10 per case
60cm (23.6"), 0.55 mL	5 French	ENSI6050	10 per case
90cm (35.4"), 0.75 mL	5 French	ENSI9050	10 per case
40cm (15.7"), 0.60 mL	6.5 French	ENSI4065	10 per case
60cm (23.6"), 0.80 mL	6.5 French	ENSI6065	10 per case
90cm (35.4"), 1.20 mL	6.5 French	ENSI9065	10 per case
40cm (15.7"), 0.90 mL	8 French	ENSI4080	10 per case
60cm (23.6"), 1.30 mL	8 French	ENSI6080	10 per case
90cm (35.4"), 1.85 mL	8 French	ENSI9080	10 per case

## Polyurethane Enteral Feeding Tubes

(For use up to 30 days)



Length/Priming Vol.	Size	Item#	Quantity
16", 0.16 mL	4 French	ENPU1640	10 per case
20", 0.18 mL	4 French	ENPU2040	10 per case
16", 0.30 mL	5 French	ENPU1650	10 per case
20", 0.35 mL	5 French	ENPU2050	10 per case
36", 0.65 mL	5 French	ENPU3650	10 per case
16", 0.50 mL	6.5 French	ENPU1665	10 per case
20", 0.60 mL	6.5 French	ENPU2065	10 per case
36", 1.10 mL	6.5 French	ENPU3665	10 per case
16", 0.80 mL	8 French	ENPU1680	10 per case
20", 1.00 mL	8 French	ENPU2080	10 per case
42", 2.00 mL	8 French	ENPU4280	10 per case
42", 3.00 mL	10 French	ENPU42100	10 per case

## Enteral Straw (for feeding preparation)


TwistLok® Oral Female Connector, Orange-Striped Standard Bore Tubing




Length/Priming Vol.	Item#	Quantity
5", 0.40 mL	TL001	100 per case


## G-Tube Extension Set (For Gastric Feeding)


TwistLok Oral Female Connector, Pinch Clamp, Enteral Feeding Direction Tag, Orange Universal Enteral Connector





Length/Priming Vol.	Item#	Quantity
12", 1.20 mL	TL002	50 per case
12", 1.60 mL	TL002Y (med. port)	50 per case

Cap	Description	Item #	Quantity
	TwistLok® Non-Vented Cap	TLCAP	500 (50 pks of 10) per case

Syringe Lid	Description	Item #	Quantity
	GraviFeed 10 mL Lid	LID10	50 per case
	GraviFeed 30 mL Lid	LID30	50 per case
	GraviFeed 60 mL Lid	LID60	50 per case

GraviFeed® Syringe	Description	Item #	Quantity
	GraviFeed 10 mL Syringe	NGF10	100 per case
	GraviFeed 30 mL Syringe	NGF30	50 per case
	GraviFeed 60 mL Syringe	NGF60	50 per case

Syringe Lid	Description	Item #	Quantity
	10 mL Plunger	PLG10	25 per case
	30 mL Plunger	PLG30	25 per case
	60 mL Plunger	PLG60	25 per case

GraviFeed Primer	Description	Item#	Quantity
	GraviFeed Primer	GF002	10 per case
	Large GraviFeed Primer	GF002L	10 per case

All extension sets, G-tube adapters, enteral straws and PVC feeding tubes are not made with DEHP or natural rubber latex.



Transforming lives, drop by drop.

Medela , 1101 Corporate Drive, McHenry, IL 60050  
Phone 888-MEDELA-8 or (815) 363-1166  
Email [accountsupport@medela.com](mailto:accountsupport@medela.com) [www.medela.com](http://www.medela.com)

Content appearance may vary from pictures.

Medela wordmark and logo, TwistLok, Every Drop Counts, Guardian Warmer, and GraviFeed are registered in the U.S. Patent and Trademark Office.

1548984 B 0418 © 2018 Medela 4550107

**medela**   
Every Drop Counts®