

Low Dose Dual Cap

Meets ISO standard for syringe accuracy when used as indicated.


The Medela Low Dose Dual Cap provides a safe, accurate, and convenient way for your pharmacy to continue the use of ISMP recommended slip-tip syringes, while allowing the entire hospital to benefit from transitioning to ENFit® enteral connectors. Medela’s Low Dose Dual Cap does not require twisting to apply the cap to the syringe. This feature results in no change to current filling practices, which may result in saved time and labor.

Benefits:

- Pharmacy can use the same oral slip-tip syringes and filling methods.
- 2-piece cap facilitates a clean connection.
- Removal of the first cap facilitates connection to ENFit feeding tubes at time of administration.
- Designed with Low Dose technology, which provides accurate dosing through the feeding tube.
- Low Dose Dual Cap is only needed when administering medication through an ENFit feeding tube.



Just 3 easy steps....



Description	General Guidance	Item #	Qty/Case (each)
Low Dose Dual Cap (LDDC) with Cap, Sterile	NICU/PICU unit	101033098	10 per pack / 50 packs per case
Low Dose Dual Cap (LDDC) with Cap, Non-Sterile	Pharmacy	101033099	500 per case (bulk)
Low Dose Dual Cap (LDDC) without Cap, Sterile	All other units	101033100	50 per case (individually wrapped)

To order or for more information, contact your Medela Sales Consultant.