

Medela DCS Disposable Fluid Collection System

SAFETY MADE EASY

Peace of mind...

- Safe, hygienic fluid collection and disposal
- Quiet operation < 50 dB(A)
- Easy handling, single-handed operation

...through smart design

- Embedded high-flow bacteria and overflow filter
- Highly burst-resistant, double-layered liners
- Automatic vacuum connection – no extra tubing



SAFETY MADE EASY

Medela's Disposable Fluid Collection System supports the safe collection and disposal of suctioned fluids.

All liners are equipped with multiple safety features such as an overflow and bacteria filter, a smoke pre-filter and a non-return valve to prevent back-flow. The hard-top lid seals precisely and immediately enable operation at low noise levels. Because the vacuum is automatically attached to the liner through the jar, no extra vacuum tubing needs to be connected with every liner exchange.

Setups with switch-over valve or serial connection enable handling of large volumes. The liners are made of non-toxic materials and can be disposed of by incineration in an environmentally friendly way.



Key features at a glance



Low noise emission significantly reduces disturbance in the operating room

- With 50 dB(A) measured at 1m distance, Medela's liners are up to four times quieter than comparable models¹.



EASY TO USE

Easy, smooth, single-handed operation

- The vacuum is immediately applied through the hard-top lids, which snap onto the jar easily and are a perfect fit.
- The large liner handles, V-shaped jars and the jar's locking clasp enable single-handed operation.
- Thanks to the automatic vacuum connection through the jar, no extra tubing needs to be connected every time



Safety by design

- Closed system with high-flow bacteria (99.999% retention¹) and overflow filter, smoke pre-filter and patient backflow valve.
- Double layered, highly burst-resistant¹ PA/PE foil ensures optimum protection.
- Optional solidifier models eliminate the risk of spilling and may reduce cost of disposal².



Full, compatible accessories range

Medela offers a full range of compatible accessories such as reusable jars, vacuum tubing, filters, specimen cups, patient tubings and carousel trolley.



References

¹ Performance data on file at Medela AG

² Depending on local disposal regulations

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Medical Vacuum Technology
for Healthcare Professionals

Please contact us or your local Medela representative for details.

 Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

Local contact